

Package ‘AffyCompatible’

September 23, 2012

Type Package

Title Affymetrix GeneChip software compatibility

Version 1.16.0

Author Martin Morgan, Robert Gentleman

Maintainer Martin Morgan <mtmorgan@fhcrc.org>

Depends R (>= 2.7.0), XML (>= 2.8-1), RCurl (>= 0.8-1), methods

Imports Biostrings

Description This package provides an interface to Affymetrix chip annotation and sample attribute files. The package allows an easy way for users to download and manage local data bases of Affymetrix NetAffx annotation files. The package also provides access to GeneChip Operating System (GCOS) and GeneChip Command Console (AGCC)-compatible sample annotation files.

License Artistic-2.0

biocViews Infrastructure, Microarray, OneChannel

R topics documented:

AffyCompatible-package	2
accessors	2
AffyCompatible-class	3
Internal-class	4
NetAffxResource	4
NetAffxResource-class	5
readArr	8
readMage	8
readNetAffx	9
readXml	10
xclass	11
xmlValue.XMLAttributeValue	12
zzz-ARRCompatibility-class	12
zzz-DTTCCompatibility-class	13
zzz-NetAffxCompatibility-class	14

Index	16
--------------	-----------

AffyCompatible-package

Affymetrix chip annotation and sample attribute management

Description

This package manages Affymetrix NetAffx annotation files, and parses MAGE and ARR sample attribute files into R data structures.

NetAffx provides a way to manage the diverse annotation files provided by Affymetrix for their arrays.

ARR files provide attribute information about samples used in an experiment. These files are created, for instance, by the GeneChip Command Console or compliant applications, or

MAGE files also provide information about samples. Support in this package is for MAGE produced by Affymetrix GCOS compliant applications, and emphasizes sample description and processing.

Details

Package details are available in the DESCRIPTION file of the package, or with the R command `packageDescription("AffyCompatible")`

Overviews of how to use the functionality of this package are available in the package vignettes. See `vignette(package="AffyCompatible")`.

Author(s)

Written by Martin Morgan.

Maintainer: Martin Morgan <martin.t.morgan@gmail.com>

accessors

Access NetAffxResource-class content

Description

These functions exist to access the information represented in objects of `linkS4class{NetAffxResource}` and related classes.

Usage

```
annotationsFile(x)
directory(x)
password(x)
user(x)
```

Arguments

x An object for which the accessor is defined.

Value

See, e.g., `help('directory,NetAffxResource-method')` for values returned by specific methods.

Author(s)

Martin Morgan

See Also

[NetAffxResource](#)

AffyCompatible-class *Virtual classes defining the "AffyCompatibility" hierarchy*

Description

The classes `AffyCompatible`, `DTTCompatibility`, `ARRCompatibility`, and `NetAffxCompatibility` define the overall class hierarchy in this package. Information about these classes is not usually relevant to the user.

Objects from the Class

These are virtual classes: no objects may be created from them.

Methods

show signature(object = "AffyCompatible"): for each slot, provide the slot name, and then invoke `show` for the slot content.

Author(s)

Martin Morgan

Examples

```
showClass(getClass("NetAffxCompatibility",
  where=getNamespace("AffyCompatible")))
```

Internal-class *(Internal) Represent 'typed associations' from XML documents*

Description

This internal class represents 'typed associates' (vectors of identically-typed objects) that appear as vectors to the R user.

Objects from the Class

Objects are created internally, during XML document parsing.

Slots

.Data: Object of class "list" containing the association objects
 .extends: Object of class "character" and length 1 describing the class each element in .Data must extend

Extends

Class "list", from data part. Class "vector", by class "list", distance 2.

Methods

[signature(x = ".TypedAssociation", i = "numeric", j = "missing"): Subset the typed association, returning a typed association with elements i.
 [signature(x = ".TypedAssociation", i = "ANY", j = "ANY"): This method is present to stop access to subsetting methods for class list; it has no user-visible functionality.
show signature(object = ".TypedAssociation"): Compactly display the typed associate, e.g., for .extends=Foo and length 5, display Foo(5).

Author(s)

Martin Morgan

NetAffxResource *Create and access a NetAffxResource database*

Description

Use this function to create and access a NetAffxResource database.

Usage

```
NetAffxResource(user = character(0),
  password = character(0),
  affxLicence = "FHCR0607",
  directory = tempdir(), ...)
```

Arguments

user	Length 1 character vector containing the Affymetrix user name authorized to access the NetAffx resource.
password	Length 1 character vector containing the Affymetrix password of the user.
affxLicence	Length 1 character vector containing the Affymetrix licence (typically application specific, i.e., the specified licence covers AffyCompatible).
directory	(Optional) length 1 character vector providing the full path to the local directory where the NetAffx database and files will be managed. Defaults to a session-specific temporary directory.
...	Additional arguments, not currently used.

Details

NetAffxResource creates a [NetAffxResource](#) object that manages a data base of NetAffx resources. The data base is located on the user file system as indicated by the directory argument. It is populated when the user requests specific files.

Users should provide a directory argument, so that the data base can be reused across sessions. This means that large files are only downloaded once.

Additional details are provided in the vignette accompanying this package.

Value

[NetAffxResource](#).

Author(s)

Martin Morgan

See Also

[NetAffxResource](#). The vignette ‘Annotations with NetAffx’.

NetAffxResource-class *Manage NetAffx annotations*

Description

This class is for accessing and retrieving NetAffx annotation files for Affymetrix chips. It requires that the user has established an account with Affymetrix. Detailed usage is described in the ‘Annotations with NetAffx’ vignette.

Objects from the Class

Objects are created by calls of the form `NetAffxResource(user="user@research.org", password="terces")`, as illustrated below and in the ‘Annotations with NetAffx’ vignette.

Slots

user: Object of class "character" containing the full user name (typically an email address) obtained when the user registers with Affymetrix.

password: Object of class "character" containing the password associated with the Affymetrix user name. The password is entered and stored in plain text, and should not be considered secure within Bioconductor.

affxLicence: Object of class "character" containing the licence, obtained from Affymetrix, for retrieving NetAffx annotations. The licence is usually software specific (e.g., for the Affy-Compatible package of Bioconductor). Users will not usually change the default value.

directory: Object of class "character", containing the local directory path where the NetAffx data base and annotation files will reside.

annotationsFile: Object of class "character" containing the file name of the NetAffx data base. Users will not usually change the default value.

affxUrl: Object of class "character" containing the URL used to download the NetAffx annotations file. Users will not usually change the default value.

Extends

Class "[NetAffxCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "NetAffxCompatibility", distance 2.

Methods

[[signature(x = "NetAffxResource", i = "numeric", j = "missing")

[[signature(x = "NetAffxResource", i = "character", j = "missing"): Retrieve an object representing all NetAffx annotations available for an array. The argument i may be a numeric index or character string corresponding to an array name returned by names(x).

[[signature(x = "NetAffxResource", i = "numeric", j = "numeric")

[[signature(x = "NetAffxResource", i = "character", j = "numeric")

[[signature(x = "NetAffxResource", i = "numeric", j = "character")

[[signature(x = "NetAffxResource", i = "character", j = "character"): Retrieve an object representing a specific annotation for an array. The argument i may be a numeric index or character string corresponding to an array name returned by names(x). The argument j may be a numeric index or character string corresponding to an annotation returned by description(x[[i]]).

affxLicence signature(x = "NetAffxResource"): the character string representation of the licence used to access the NetAffx site.

affxUrl signature(x = "NetAffxResource"): the url used to access the NetAffx site.

annotationsFile signature(x = "NetAffxResource"): the file name of NetAffx data base, as stored on the local disk of the user.

directory signature(x = "NetAffxResource"): the directory in which the NetAffx data base and annotation files are stored on the local disk of the user.

names signature(x = "NetAffxResource"): the names of the arrays for which annotations are available.

password signature(x = "NetAffxResource"): The password associated with the user name obtained when the user registered with Affymetrix. The password is entered and stored as plain text, and is not secure.

readAnnotation signature(netAffxResource = "NetAffxResource", array = "numeric", annotation = "n

readAnnotation signature(netAffxResource = "NetAffxResource", array = "numeric", annotation = "c

readAnnotation signature(netAffxResource = "NetAffxResource", array = "character", annotation =

readAnnotation signature(netAffxResource = "NetAffxResource", array = "character", annotation =

readAnnotation signature(netAffxResource = "NetAffxResource", array = "missing", annotation = "A

Read a specific annotation file as an R object. For the first four methods, array, if present, indexes or names the array while annotation indexes or names the annotation. All methods also accept content = "logical" to indicate whether the content of the file should be returned if possible (content=TRUE, default) or only the path to the local file (content=FALSE) and update = "logical" indicating whether the file should be read from disk if available (update=FALSE, default) or retrieved from the NetAffx web site (update=TRUE). readAnnotation parses CSV, Tabular, and PSI files to data frames, and reads FASTA files to a list of description / sequence pairs (see Biostrings::readFASTA).

The ... argument is passed to corresponding read functions. CSV files are parsed using read.csv(conn, ...), Tabular files with read.delim(conn, ...), PSI files with read.delim(conn, header=) and FASTA with readFASTA(conn, ...).

show signature(object = "NetAffxResource"): compactly show the NetAffxResource object.

user signature(x = "NetAffxResource"): The user name (typically email address) obtained when the user registers with the Affymetrix web site.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Annotations with NetAffx'. Affymetrix user login site <https://www.affymetrix.com/site/login/login.affx>.

.Rprofile for ways to automatically set variables (such as password, and directory) for each session.

Examples

```
## password must be obtained by the user
rsrc <- NetAffxResource(user="mtmorgan@fhcrc.org", password=password)
head(names(rsrc))
affxDescription(rsrc[["Bovine"]])
## paradigm to display all elements of 'annotation'
sapply(affxAnnotation(rsrc[["Bovine"]]), force)
## obtain and read a single annotation
anno <- rsrc[["Porcine", "BLASTP annotations, CSV Format"]]
df <- readAnnotation(rsrc, annotation=anno)
```

readArr

Read ARR sample attribute files

Description

This functions reads the content of Affymetrix ARR sample attribute files created by GeneChip(R) Command Console. Objects are represented as R classes. Navigation and extraction of data are described in the vignette ‘Retrieving MAGE and ARR sample attributes’.

Usage

```
readArr(fls, ...)
```

Arguments

<code>fls</code>	A character vector containing the paths to files to be read.
<code>...</code>	Additional arguments, currently unused.

Value

An object or list of objects of class `ArrayFileSet`, depending on whether one or several file paths are provided.

Author(s)

Martin Morgan

See Also

The vignette ‘Retrieving MAGE and ARR sample attributes’ for detailed illustration of use. [xclass](#) and [readXml](#) for advanced access methods.

Examples

```
dir <- system.file("extdata", "ARR", package="AffyCompatible")
fls <- list.files(dir, full=TRUE)
readArr(fls[[1]])
readArr(fls)
```

readMage*Read DTT (MAGE) sample attribute files*

Description

This functions reads the content of Affymetrix DTT (MAGE) sample attribute files created by GCOS. Objects are represented as R classes. Navigation and extraction of data are described in the vignette ‘Retrieving MAGE and ARR sample attributes’.

Usage

```
readMage(fls, ...)
```

Arguments

`fls` A character vector containing the paths to files to be read.
`...` Additional arguments, currently unused.

Value

An object or list of objects of class `ArrayFileSet`, depending on whether one or several file paths are provided.

Author(s)

Martin Morgan

See Also

The vignette ‘Retrieving MAGE and ARR sample attributes’ for detailed illustration of use. [xclass](#) and [readXml](#) for advanced access methods.

Examples

```
dir <- system.file("extdata", "DTT", package="AffyCompatible")
fls <- list.files(dir, full=TRUE)
basename(fls)
readArr(fls[[1]])
```

readNetAffx

Update or read a NetAffx data base

Description

This function provides a way to read or update a NetAffx data base. The usual use for this is to force the data base description to be updated (by accessing the Affymetrix web site) or to retrieve the entire list of NetAffx resource data base. A more normal use case is to

Usage

```
readNetAffx(netAffxResource, update = FALSE, ...)
```

Arguments

`netAffxResource` An object of class [NetAffxResource](#) containing the path to the local data base, and permissions for access the Affymetrix NetAffx resource.
`update` An object of `logical(1)` indicating whether the data base is to be updated (TRUE) or not.
`...` Additional arguments, not currently used.

Value

An object representing the annotations found; see [NetAffxResource](#) for additional details.

Author(s)

Martin Morgan

Examples

```
## retrieve all annotations
## Not run:
example("NetAffxResource-class")
readNetAffx(rsrc)

## End(Not run)
```

readXml	<i>(Advanced) Read XML files of Affymetrix sample attributes or annotation</i>
---------	--

Description

This functions reads the content of Affymetrix DTT (MAGE) or ARR sample files, or AffxNetAffx-AnnotFile array annotation files, return an internal representation of the underlying XML. The return value can be used with [xclass](#) to rapidly create R objects representing deeply nested nodes of the document. The return value can also be used with `xpathApply` in the XML package for access to character-based representation of the data.

Usage

```
readXml(f1, ...)
```

Arguments

f1	A length-1 character vector giving the full path to the file to be parsed.
...	Additional arguments, passed to <code>xmlTreeParse</code> .

Value

An object of class `XMLInternalDocument`.

Author(s)

Martin Morgan

See Also

The vignettes ‘Retrieving MAGE and ARR sample attributes’ and ‘Annotation retrieval with `NetAffxResource`’ for detailed illustration of use.

[xclass](#) for instantiating objects from `xpath` queries.

Examples

```
mageDir <- system.file("extdata", "DTT", package="AffyCompatible")
mageFiles <- list.files(mageDir, full=TRUE)
xml <- readXml(mageFiles[[1]])
xclass(xml, "//MeasuredBioAssay")[[1]]
```

xclass

(Advanced) Construct R objects from XML document queries

Description

This function takes an xml object, created from [readXml](#), and an xpath query referencing an element of the document, and returns an R object representing the element.

The xpath query must retrieve an element, not an attribute.

Usage

```
xclass(doc, xpathq, prefix = character(0))
```

Arguments

doc	An XMLInternalDocument retrieved using readXml .
xpathq	A length-1 character string representing an xpath query into the document.
prefix	A prefix used to map, when necessary, between the XML class name and the corresponding R class. For NetAffx-based classes, the prefix is 'Affx'; no prefix is required for DTT or ARR documents.

Value

An R object corresponding to the class(es) referenced by the xpath query.

Author(s)

Martin Morgan

References

<http://www.w3.org/TR/xpath> provides a very useful description of xpath. Additional references are in the vignette 'Retrieving MAGE and ARR sample attributes'.

See Also

The vignettes 'Retrieving MAGE and ARR sample attributes' and 'Annotation retrieval with NetAffxResource' for detailed illustration of use.

[readXml](#) for appropriately reading XML documents.

Examples

```
example(readXml)
```

xmlValue.XMLAttributeValue

Additional functions defined in AffyCompatible

Description

This page documents additional, miscellaneous, functions defined by AffyCompatible

Usage

```
xmlValue.XMLAttributeValue(x, ignoreComments = FALSE, recursive = TRUE, encoding=XML:::CE_NATIVE)
```

Arguments

x	An object of class XMLAttributeValue.
ignoreComments	Logical; ignored by the xmlValue.XMLAttributeValue method
recursive	ignored
encoding	File encoding

Value

A character vector (usually length 1) containing the character representation of the XML attribute.

zzz-ARRCompatibility-class

(Auto-generated) classes and methods for navigating ARR sample attributes.

Description

These classes and methods facilitate navigation of the Affymetrix ARR sample attribute files. ARR files are generated by the Affymetrix AGCC (GeneChip(R) Command Console) application. Specifically, the classes are R representations of the Affymetrix classes defined in ArraySetAndTemplate-File.dtd. The 'show' methods for these class indicate how objects are to be accessed; details are in the 'Retrieving MAGE and ARR sample attributes' vignette.

Objects from the Class

Objects are created internally, by parsing XML retrieved from ARR, and are not created directly by the user.

Slots

Objects have a variety of slots. To discover slot names and content type, see the slotNames portion of the example, below.

Extends

Class "[ARRCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "ARRCompatibility", distance 2.

Methods

'Accessor' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form accessor: content. To retrieve the content of an instance object, invoke accessor(object).

Some objects contains several instances of another object. These 'typed associations' are represented by the [.TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using [, [[]]. A useful paradigm for displaying all instances is illustrated in the vignette.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Retrieving MAGE and ARR sample attributes'.

Examples

```
## defined and auto-generated classes
names(slot(getClass("ARRCompatibility",
                  where=getNamespace("AffyCompatible")),
       "subclasses"))
```

zzz-DTTCompatibility-class

(Auto-generated) classes and methods for navigating DTT sample attributes.

Description

These classes and methods facilitate navigation of the Affymetrix DTT sample attribute files. DTT files are generated by the Affymetrix GCOS application. Specifically, the classes are R representations of the Affymetrix classes defined in MAGE-ML.dtd. The 'show' methods for these class indicate how objects are to be accessed; details are in the 'Retrieving MAGE and ARR sample attributes' vignette.

Objects from the Class

Objects are created internally, by parsing XML retrieved from ARR, and are not created directly by the user.

Slots

Objects have a variety of slots. To discover slot names and content type, see the slotNames portion of the example, below.

Extends

Class "[DTTCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "[DTTCompatibility](#)", distance 2.

Methods

'Accessor' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form accessor: content. To retrieve the content of an instance object, invoke accessor(object).

Some objects contains several instances of another object. These 'typed associations' are represented by the `.TypedAssociation` class, and are navigated or subset as one would a standard R vector, e.g., using `[`, `[[`. A useful paradigm for displaying all instances is illustrated in the vignette.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Retrieving MAGE and ARR sample attributes'.

Examples

```
## defined and auto-generated classes
names(slot(getClass("DTTCompatibility",
                  where=getNamespace("AffyCompatible")),
       "subclasses"))
```

zzz-NetAffxCompatibility-class

(Auto-generated) classes and methods for navigating NetAffx resources

Description

These classes and methods facilitate navigation of the NetAffx resource data base. They are R representations of the Affymetrix classes defined in NetAffxAnnotFileList.dtd. The 'show' methods for these class indicate how objects are to be accessed; details are in the 'Annotations with NetAffx' vignette.

The auto-generated classes are discoverable as indicated in the example below.

Objects from the Class

Objects are created internally, by parsing XML retrieved from the local NetAffx data base or from the Affymetrix web site.

Slots

Objects have a variety of slots. To discover slot names and content type, see the slotNames portion of the example, below.

Extends

Class "`NetAffxCompatibility`", directly. Class "`AffyCompatible`", by class "`NetAffxCompatibility`", distance 2.

Methods

'Read' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form accessor: content. To retrieve the content of an instance object, invoke accessor(object).

Some objects contains several instances of another object. These 'typed associations' are represented by the [.TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using [, [[]. A useful paradigm for displaying all instances is illustrated on the [NetAffxResource](#) page.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Annotations with NetAffx'. Affymetrix user login site <https://www.affymetrix.com/site/login/login.affx>.

[NetAffxResource](#) for creating and managing the NetAffx data base.

[NetAffxResource](#) for navigating the data base.

Examples

```
## defined and auto-generated classes
getClass("NetAffxCompatibility", where=getNamespace("AffyCompatible"))

## slots of an auto-generated class
getSlots(getClass("AffxNetAffxAnnotFileList",
  where=getNamespace("AffyCompatible")))

## Not run:
example("NetAffxrResource-class")

## End(Not run)
```

Index

- *Topic **classes**
 - AffyCompatible-class, 3
 - Internal-class, 4
 - NetAffxResource-class, 5
 - readArr, 8
 - readMage, 8
 - zzz-ARRCompatibility-class, 12
 - zzz-DTTCompatibility-class, 13
 - zzz-NetAffxCompatibility-class, 14
- *Topic **database**
 - NetAffxResource, 4
 - readNetAffx, 9
- *Topic **manip**
 - readXml, 10
 - xclass, 11
 - xmlValue.XMLAttributeValue, 12
- *Topic **methods**
 - accessors, 2
- *Topic **package**
 - AffyCompatible-package, 2
 - .Rprofile, 7
 - .TypedAssociation, 13–15
 - .TypedAssociation-class
 - (Internal-class), 4
 - [, .TypedAssociation, ANY, ANY-method
 - (Internal-class), 4
 - [, .TypedAssociation, numeric, missing, ANY-method
 - (Internal-class), 4
 - [, .TypedAssociation, numeric, missing-method
 - (Internal-class), 4
 - [[, NetAffxResource, character, character-method
 - (NetAffxResource-class), 5
 - [[, NetAffxResource, character, missing-method
 - (NetAffxResource-class), 5
 - [[, NetAffxResource, character, numeric-method
 - (NetAffxResource-class), 5
 - [[, NetAffxResource, numeric, character-method
 - (NetAffxResource-class), 5
 - [[, NetAffxResource, numeric, missing-method
 - (NetAffxResource-class), 5
 - [[, NetAffxResource, numeric, numeric-method
 - (NetAffxResource-class), 5
- accession (zzz-DTTCompatibility-class), 13
- accession, DatabaseEntry-method
 - (zzz-DTTCompatibility-class), 13
- accessions_assnlist
 - (zzz-DTTCompatibility-class), 13
- accessions_assnlist, BibliographicReference-method
 - (zzz-DTTCompatibility-class), 13
- Accessions_assnlist-class
 - (zzz-DTTCompatibility-class), 13
- accessionVersion
 - (zzz-DTTCompatibility-class), 13
- accessionVersion, DatabaseEntry-method
 - (zzz-DTTCompatibility-class), 13
- accessors, 2
- action (zzz-DTTCompatibility-class), 13
- action, Audit-method
 - (zzz-DTTCompatibility-class), 13
- action_assn
 - (zzz-DTTCompatibility-class), 13
- action_assn, Treatment-method
 - (zzz-DTTCompatibility-class), 13
- Action_assn-class
 - (zzz-DTTCompatibility-class), 13
- actionMeasurement_assn
 - (zzz-DTTCompatibility-class), 13
- actionMeasurement_assn, Treatment-method
 - (zzz-DTTCompatibility-class), 13
- ActionMeasurement_assn-class
 - (zzz-DTTCompatibility-class), 13
- activityDate

- (zzz-DTTCompatibility-class),
13
- activityDate, ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- address (zzz-DTTCompatibility-class), 13
- address, Organization-method
(zzz-DTTCompatibility-class),
13
- address, Person-method
(zzz-DTTCompatibility-class),
13
- adjustments_assnlist
(zzz-DTTCompatibility-class),
13
- adjustments_assnlist, ArrayManufactureDeviation-method
(zzz-DTTCompatibility-class),
13
- Adjustments_assnlist-class
(zzz-DTTCompatibility-class),
13
- affiliation_assnref
(zzz-DTTCompatibility-class),
13
- affiliation_assnref, Person-method
(zzz-DTTCompatibility-class),
13
- Affiliation_assnref-class
(zzz-DTTCompatibility-class),
13
- affxAnnotation
(zzz-NetAffxCompatibility-class),
14
- affxAnnotation, AffxArray-method
(zzz-NetAffxCompatibility-class),
14
- AffxAnnotation-class
(zzz-NetAffxCompatibility-class),
14
- affxArray
(zzz-NetAffxCompatibility-class),
14
- affxArray, AffxNetAffxAnnotFileList-method
(zzz-NetAffxCompatibility-class),
14
- affxArray, Array_assnlist-method
(zzz-DTTCompatibility-class),
13
- AffxArray-class
(zzz-NetAffxCompatibility-class),
14
- affxCategory
(zzz-DTTCompatibility-class),
13
- affxCategory, OntologyEntry-method
(zzz-DTTCompatibility-class),
13
- affxComment
(zzz-NetAffxCompatibility-class),
14
- affxComment, AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
14
- affxComment, PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- AffxComment-class
(zzz-NetAffxCompatibility-class),
14
- affxCompression
(zzz-NetAffxCompatibility-class),
14
- affxCompression, AffxUrl-method
(zzz-NetAffxCompatibility-class),
14
- affxCrc32
(zzz-NetAffxCompatibility-class),
14
- affxCrc32, AffxFile-method
(zzz-NetAffxCompatibility-class),
14
- affxCrc32, AffxUrl-method
(zzz-NetAffxCompatibility-class),
14
- affxDate
(zzz-NetAffxCompatibility-class),
14
- affxDate, AffxFile-method
(zzz-NetAffxCompatibility-class),
14
- affxDate, Audit-method
(zzz-DTTCompatibility-class),
13
- affxDescription
(zzz-NetAffxCompatibility-class),
14
- affxDescription, AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
14
- affxDescription, AffxArray-method
(zzz-NetAffxCompatibility-class),
14
- affxEnd (zzz-DTTCompatibility-class), 13
- affxEnd, CompositePosition-method

- (zzz-DTTCompatibility-class),
13
- affxEnd, ReporterPosition-method
(zzz-DTTCompatibility-class),
13
- affxEnd, SequencePosition-method
(zzz-DTTCompatibility-class),
13
- affxFile
(zzz-NetAffxCompatibility-class),
14
- affxFile, AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
14
- AffxFile-class
(zzz-NetAffxCompatibility-class),
14
- affxImage (zzz-DTTCompatibility-class),
13
- affxImage, PhysicalBioAssayData_assnlist-method
(zzz-DTTCompatibility-class),
13
- affxLicence
(zzz-NetAffxCompatibility-class),
14
- affxLicence, NetAffxResource-method
(NetAffxResource-class), 5
- affxName
(zzz-NetAffxCompatibility-class),
14
- affxName, AffxArray-method
(zzz-NetAffxCompatibility-class),
14
- affxName, AffxFile-method
(zzz-NetAffxCompatibility-class),
14
- AffxNetAffxAnnotFileList-class
(zzz-NetAffxCompatibility-class),
14
- affxOrder (zzz-DTTCompatibility-class),
13
- affxOrder, BioDataCube-method
(zzz-DTTCompatibility-class),
13
- affxOrder, Treatment-method
(zzz-DTTCompatibility-class),
13
- affxPerson
(zzz-DTTCompatibility-class),
13
- affxPerson, Contact_assnlist-method
(zzz-DTTCompatibility-class),
13
- affxPerson, Person_assnlist-method
(zzz-DTTCompatibility-class),
13
- affxRow (zzz-DTTCompatibility-class), 13
- affxRow, FeatureLocation-method
(zzz-DTTCompatibility-class),
13
- affxRow, Zone-method
(zzz-DTTCompatibility-class),
13
- affxSequence
(zzz-DTTCompatibility-class),
13
- affxSequence, BioSequence-method
(zzz-DTTCompatibility-class),
13
- affxSize
(zzz-NetAffxCompatibility-class),
14
- affxSize, AffxFile-method
(zzz-NetAffxCompatibility-class),
14
- affxSize, AffxUrl-method
(zzz-NetAffxCompatibility-class),
14
- affxStart (zzz-DTTCompatibility-class),
13
- affxStart, CompositePosition-method
(zzz-DTTCompatibility-class),
13
- affxStart, ReporterPosition-method
(zzz-DTTCompatibility-class),
13
- affxStart, SequencePosition-method
(zzz-DTTCompatibility-class),
13
- affxStatus
(zzz-NetAffxCompatibility-class),
14
- affxStatus, AffxNetAffxAnnotFileList-method
(zzz-NetAffxCompatibility-class),
14
- AffxStatus-class
(zzz-NetAffxCompatibility-class),
14
- affxText (zzz-DTTCompatibility-class),
13
- affxText, Description-method
(zzz-DTTCompatibility-class),
13
- affxText, Protocol-method

- (zzz-DTTCompatibility-class),
13
- affxTitle (zzz-DTTCompatibility-class),
13
- affxTitle, BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- affxTitle, Protocol-method
(zzz-DTTCompatibility-class),
13
- affxType
(zzz-NetAffxCompatibility-class),
14
- affxType, AffxAnnotation-method
(zzz-NetAffxCompatibility-class),
14
- affxType, AffxFile-method
(zzz-NetAffxCompatibility-class),
14
- affxUrl
(zzz-NetAffxCompatibility-class),
14
- affxUrl, AffxFile-method
(zzz-NetAffxCompatibility-class),
14
- affxUrl, NetAffxResource-method
(NetAffxResource-class), 5
- AffxUrl-class
(zzz-NetAffxCompatibility-class),
14
- affyBarcode
(zzz-ARRCompatibility-class),
12
- affyBarcode, PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- AffyCompatible, 6, 12–14
- AffyCompatible
(AffyCompatible-package), 2
- AffyCompatible-class, 3
- AffyCompatible-package, 2
- analysisResults_assnreflist
(zzz-DTTCompatibility-class),
13
- analysisResults_assnreflist, Experiment-method
(zzz-DTTCompatibility-class),
13
- AnalysisResults_assnreflist-class
(zzz-DTTCompatibility-class),
13
- annotations_assnlist
(zzz-DTTCompatibility-class),
13
- annotations_assnlist, Description-method
(zzz-DTTCompatibility-class),
13
- annotations_assnlist, ExperimentalFactor-method
(zzz-DTTCompatibility-class),
13
- Annotations_assnlist-class
(zzz-DTTCompatibility-class),
13
- annotationsFile (accessors), 2
- annotationsFile, NetAffxResource-method
(NetAffxResource-class), 5
- Array-class
(zzz-DTTCompatibility-class),
13
- array_assnlist
(zzz-DTTCompatibility-class),
13
- array_assnlist, Array_package-method
(zzz-DTTCompatibility-class),
13
- Array_assnlist-class
(zzz-DTTCompatibility-class),
13
- array_assnref
(zzz-DTTCompatibility-class),
13
- array_assnref, BioAssayCreation-method
(zzz-DTTCompatibility-class),
13
- array_assnref, Hybridization-method
(zzz-DTTCompatibility-class),
13
- Array_assnref-class
(zzz-DTTCompatibility-class),
13
- array_package
(zzz-DTTCompatibility-class),
13
- array_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- Array_package-class
(zzz-DTTCompatibility-class),
13
- array_ref (zzz-DTTCompatibility-class),
13
- array_ref, Array_assnref-method
(zzz-DTTCompatibility-class),
13
- array_ref, Arrays_assnreflist-method

- (zzz-DTTCompatibility-class),
13
- Array_ref-class
(zzz-DTTCompatibility-class),
13
- arrayAttribute
(zzz-ARRCompatibility-class),
12
- arrayAttribute,PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- ArrayAttribute-class
(zzz-ARRCompatibility-class),
12
- arrayDesign
(zzz-DTTCompatibility-class),
13
- arrayDesign,ArrayDesign_assnlist-method
(zzz-DTTCompatibility-class),
13
- ArrayDesign-class
(zzz-DTTCompatibility-class),
13
- arrayDesign_assnlist
(zzz-DTTCompatibility-class),
13
- arrayDesign_assnlist,ArrayDesign_package-method
(zzz-DTTCompatibility-class),
13
- ArrayDesign_assnlist-class
(zzz-DTTCompatibility-class),
13
- arrayDesign_assnref
(zzz-DTTCompatibility-class),
13
- arrayDesign_assnref,Array-method
(zzz-DTTCompatibility-class),
13
- ArrayDesign_assnref-class
(zzz-DTTCompatibility-class),
13
- arrayDesign_package
(zzz-DTTCompatibility-class),
13
- arrayDesign_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- ArrayDesign_package-class
(zzz-DTTCompatibility-class),
13
- arrayDesign_ref
(zzz-DTTCompatibility-class),
13
- arrayDesign_ref,ArrayDesign_assnref-method
(zzz-DTTCompatibility-class),
13
- ArrayDesign_ref-class
(zzz-DTTCompatibility-class),
13
- arrayGroup
(zzz-DTTCompatibility-class),
13
- arrayGroup,ArrayGroup_assnlist-method
(zzz-DTTCompatibility-class),
13
- ArrayGroup-class
(zzz-DTTCompatibility-class),
13
- arrayGroup_assnlist
(zzz-DTTCompatibility-class),
13
- arrayGroup_assnlist,Array_package-method
(zzz-DTTCompatibility-class),
13
- ArrayGroup_assnlist-class
(zzz-DTTCompatibility-class),
13
- arrayGroup_assnref
(zzz-DTTCompatibility-class),
13
- arrayGroup_assnref,Array-method
(zzz-DTTCompatibility-class),
13
- ArrayGroup_assnref-class
(zzz-DTTCompatibility-class),
13
- arrayGroup_ref
(zzz-DTTCompatibility-class),
13
- arrayGroup_ref,ArrayGroup_assnref-method
(zzz-DTTCompatibility-class),
13
- ArrayGroup_ref-class
(zzz-DTTCompatibility-class),
13
- arrayIdentifier
(zzz-DTTCompatibility-class),
13
- arrayIdentifier,Array-method
(zzz-DTTCompatibility-class),
13
- arrayManufacture
(zzz-DTTCompatibility-class),
13

- arrayManufacture,ArrayManufacture_assnlist-method 12
(zzz-DTTCompatibility-class), 13
- ArrayManufacture-class
(zzz-DTTCompatibility-class), 13
- arrayManufacture_assnlist
(zzz-DTTCompatibility-class), 13
- arrayManufacture_assnlist,Array_package-method
(zzz-DTTCompatibility-class), 13
- ArrayManufacture_assnlist-class
(zzz-DTTCompatibility-class), 13
- arrayManufacture_ref
(zzz-DTTCompatibility-class), 13
- arrayManufacture_ref,Information_assnref-method
(zzz-DTTCompatibility-class), 13
- ArrayManufacture_ref-class
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviation
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviation,ArrayManufactureDeviations_assnlist-method
(zzz-DTTCompatibility-class), 13
- ArrayManufactureDeviation-class
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviations_assnlist
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviations_assnlist,Array-method
(zzz-DTTCompatibility-class), 13
- ArrayManufactureDeviations_assnlist-class
(zzz-DTTCompatibility-class), 13
- arrayManufacturers_assnreflist
(zzz-DTTCompatibility-class), 13
- arrayManufacturers_assnreflist,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- ArrayManufacturers_assnreflist-class
(zzz-DTTCompatibility-class), 13
- arrayName (zzz-ARRCompatibility-class),
arrayName,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- arrays_assnreflist
(zzz-DTTCompatibility-class), 13
- arrays_assnreflist,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- arrays_assnreflist,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- Arrays_assnreflist-class
(zzz-DTTCompatibility-class), 13
- ArraySetFile-class
(zzz-ARRCompatibility-class), 12
- arraySpacingX
(zzz-DTTCompatibility-class), 13
- arraySpacingX,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- arraySpacingY
(zzz-DTTCompatibility-class), 13
- arraySpacingY,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- arrayXOrigin
(zzz-DTTCompatibility-class), 13
- arrayXOrigin,Array-method
(zzz-DTTCompatibility-class), 13
- arrayYOrigin
(zzz-DTTCompatibility-class), 13
- arrayYOrigin,Array-method
(zzz-DTTCompatibility-class), 13
- ARRCompatibility, 12
- ARRCompatibility-class
(AffyCompatible-class), 3
- associations_assnlist
(zzz-DTTCompatibility-class), 13
- associations_assnlist,OntologyEntry-method
(zzz-DTTCompatibility-class), 13

- Associations_assnlist-class
(zzz-DTTCompatibility-class),
[13](#)
- audit (zzz-DTTCompatibility-class), [13](#)
- audit, AuditTrail_assnlist-method
(zzz-DTTCompatibility-class),
[13](#)
- Audit-class
(zzz-DTTCompatibility-class),
[13](#)
- auditAndSecurity_package
(zzz-DTTCompatibility-class),
[13](#)
- auditAndSecurity_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
[13](#)
- AuditAndSecurity_package-class
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Array-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, ArrayDesign-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, ArrayGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, ArrayManufacture-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Audit-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BibliographicReference-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioAssayCreation-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioAssayDimension-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioAssayMap-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioAssayTreatment-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioSample-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioSequence-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, BioSource-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Channel-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, CompositeGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, CompositeSequence-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Compound-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Database-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, DerivedSignal-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Description-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Error-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, ExpectedValue-method
(zzz-DTTCompatibility-class),
[13](#)
- auditTrail_assnlist, Experiment-method

- (zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ExperimentalFactor-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ExperimentDesign-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,FactorValue-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Failed-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Feature-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,FeatureDimension-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Fiducial-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Hardware-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,HardwareApplication-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Hybridization-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Image-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ManufactureLIMS-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Node-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,NodeContents-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Organization-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Parameter-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Person-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Protocol-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,PValue-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
13

- auditTrail_assnlist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Ratio-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Reporter-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ReporterDimension-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Security-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,SecurityGroup-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,SeqFeature-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Software-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,SoftwareApplication-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Transformation-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Treatment-method
(zzz-DTTCompatibility-class),
13
- auditTrail_assnlist,Zone-method
(zzz-DTTCompatibility-class),
13
- AuditTrail_assnlist-class
(zzz-DTTCompatibility-class),
13
- authors (zzz-DTTCompatibility-class), 13
- authors,BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- barcode (zzz-DTTCompatibility-class), 13
- barcode,ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- basis (zzz-DTTCompatibility-class), 13
- basis,SeqFeature-method
(zzz-DTTCompatibility-class),
13
- bibliographicReference
(zzz-DTTCompatibility-class),
13
- bibliographicReference,BibliographicReferences_assnlist
(zzz-DTTCompatibility-class),
13
- BibliographicReference-class
(zzz-DTTCompatibility-class),
13
- bibliographicReferences_assnlist
(zzz-DTTCompatibility-class),
13
- bibliographicReferences_assnlist,Description-method
(zzz-DTTCompatibility-class),
13
- BibliographicReferences_assnlist-class
(zzz-DTTCompatibility-class),
13
- bioAssay_assnlist
(zzz-DTTCompatibility-class),
13
- bioAssay_assnlist,BioAssay_package-method
(zzz-DTTCompatibility-class),
13
- BioAssay_assnlist-class
(zzz-DTTCompatibility-class),
13
- bioAssay_assnref
(zzz-DTTCompatibility-class),
13
- bioAssay_assnref,BioAssayTuple-method
(zzz-DTTCompatibility-class),
13
- BioAssay_assnref-class
(zzz-DTTCompatibility-class),
13
- bioAssay_package
(zzz-DTTCompatibility-class),
13
- bioAssay_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- BioAssay_package-class
(zzz-DTTCompatibility-class),

- 13
- BioAssay_ref-class
(zzz-DTTCompatibility-class),
13
- bioAssayCreation
(zzz-DTTCompatibility-class),
13
- bioAssayCreation,BioAssayCreation_assn-method
(zzz-DTTCompatibility-class),
13
- BioAssayCreation-class
(zzz-DTTCompatibility-class),
13
- bioAssayCreation_assn
(zzz-DTTCompatibility-class),
13
- bioAssayCreation_assn,PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- BioAssayCreation_assn-class
(zzz-DTTCompatibility-class),
13
- bioAssayData_assnlist
(zzz-DTTCompatibility-class),
13
- bioAssayData_assnlist,BioAssayData_package-method
(zzz-DTTCompatibility-class),
13
- BioAssayData_assnlist-class
(zzz-DTTCompatibility-class),
13
- bioAssayData_assnreflist
(zzz-DTTCompatibility-class),
13
- bioAssayData_assnreflist,Experiment-method
(zzz-DTTCompatibility-class),
13
- BioAssayData_assnreflist-class
(zzz-DTTCompatibility-class),
13
- bioAssayData_package
(zzz-DTTCompatibility-class),
13
- bioAssayData_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- BioAssayData_package-class
(zzz-DTTCompatibility-class),
13
- BioAssayData_ref-class
(zzz-DTTCompatibility-class),
13
- bioAssayDataCluster
(zzz-DTTCompatibility-class),
13
- bioAssayDataCluster,BioAssayDataCluster_assnlist-method
(zzz-DTTCompatibility-class),
13
- BioAssayDataCluster-class
(zzz-DTTCompatibility-class),
13
- bioAssayDataCluster_assnlist
(zzz-DTTCompatibility-class),
13
- bioAssayDataCluster_assnlist,HigherLevelAnalysis_package
(zzz-DTTCompatibility-class),
13
- BioAssayDataCluster_assnlist-class
(zzz-DTTCompatibility-class),
13
- bioAssayDataCluster_ref
(zzz-DTTCompatibility-class),
13
- bioAssayDataCluster_ref,AnalysisResults_assnreflist-method
(zzz-DTTCompatibility-class),
13
- BioAssayDataCluster_ref-class
(zzz-DTTCompatibility-class),
13
- bioAssayDataSources_assnreflist
(zzz-DTTCompatibility-class),
13
- bioAssayDataSources_assnreflist,Transformation-method
(zzz-DTTCompatibility-class),
13
- BioAssayDataSources_assnreflist-class
(zzz-DTTCompatibility-class),
13
- bioAssayDimension
(zzz-DTTCompatibility-class),
13
- bioAssayDimension,BioAssayDimension_assnlist-method
(zzz-DTTCompatibility-class),
13
- BioAssayDimension-class
(zzz-DTTCompatibility-class),
13
- bioAssayDimension_assnlist
(zzz-DTTCompatibility-class),
13
- bioAssayDimension_assnlist,BioAssayData_package-method
(zzz-DTTCompatibility-class),
13
- BioAssayDimension_assnlist-class

- (zzz-DTTCompatibility-class), 13
- bioAssayDimension_assnref (zzz-DTTCompatibility-class), 13
- bioAssayDimension_assnref, DerivedBioAssayData-method (zzz-DTTCompatibility-class), 13
- bioAssayDimension_assnref, MeasuredBioAssayData-method (zzz-DTTCompatibility-class), 13
- bioAssayDimension_assnref, NodeContents-method (zzz-DTTCompatibility-class), 13
- BioAssayDimension_assnref-class (zzz-DTTCompatibility-class), 13
- bioAssayDimension_ref (zzz-DTTCompatibility-class), 13
- bioAssayDimension_ref, BioAssayDimension_assnref-method (zzz-DTTCompatibility-class), 13
- BioAssayDimension_ref-class (zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_assnreflist (zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_assnreflist, DerivedBioAssay-method (zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_assnreflist, MeasuredBioAssay-method (zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_assnreflist, PhysicalBioAssay-method (zzz-DTTCompatibility-class), 13
- BioAssayFactorValues_assnreflist-class (zzz-DTTCompatibility-class), 13
- bioAssayMap (zzz-DTTCompatibility-class), 13
- bioAssayMap, BioAssayMap_assnlist-method (zzz-DTTCompatibility-class), 13
- BioAssayMap-class (zzz-DTTCompatibility-class), 13
- bioAssayMap_assnlist (zzz-DTTCompatibility-class), 13
- bioAssayMap_assnlist, BioAssayData_package-method (zzz-DTTCompatibility-class), 13
- BioAssayMap_assnlist-class
- bioAssayMap_ref (zzz-DTTCompatibility-class), 13
- bioAssayMap_ref, BioAssayMaps_assnreflist-method (zzz-DTTCompatibility-class), 13
- bioAssayMap_ref, DerivedBioAssayMap_assnreflist-method (zzz-DTTCompatibility-class), 13
- BioAssayMap_ref-class (zzz-DTTCompatibility-class), 13
- bioAssayMapping (zzz-DTTCompatibility-class), 13
- bioAssayMapping, BioAssayMapping_assn-method (zzz-DTTCompatibility-class), 13
- BioAssayMapping-class (zzz-DTTCompatibility-class), 13
- bioAssayMapping_assn (zzz-DTTCompatibility-class), 13
- bioAssayMapping_assn, Transformation-method (zzz-DTTCompatibility-class), 13
- BioAssayMapping_assn-class (zzz-DTTCompatibility-class), 13
- bioAssayMaps_assnreflist (zzz-DTTCompatibility-class), 13
- bioAssayMaps_assnreflist, BioAssayMapping-method (zzz-DTTCompatibility-class), 13
- BioAssayMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- bioAssayMapTarget_assnref (zzz-DTTCompatibility-class), 13
- bioAssayMapTarget_assnref, BioAssayMap-method (zzz-DTTCompatibility-class), 13

- BioAssayMapTarget_assnref-class (zzz-DTTCompatibility-class), 13
- bioAssays_assnreflist (zzz-DTTCompatibility-class), 13
- bioAssays_assnreflist,BioAssayDimension-method (zzz-DTTCompatibility-class), 13
- bioAssays_assnreflist,Experiment-method (zzz-DTTCompatibility-class), 13
- BioAssays_assnreflist-class (zzz-DTTCompatibility-class), 13
- bioAssayTreatment (zzz-DTTCompatibility-class), 13
- bioAssayTreatment,BioAssayTreatments_assnlist-method (zzz-DTTCompatibility-class), 13
- BioAssayTreatment-class (zzz-DTTCompatibility-class), 13
- bioAssayTreatments_assnlist (zzz-DTTCompatibility-class), 13
- bioAssayTreatments_assnlist,PhysicalBioAssay-method (zzz-DTTCompatibility-class), 13
- BioAssayTreatments_assnlist-class (zzz-DTTCompatibility-class), 13
- bioAssayTuple (zzz-DTTCompatibility-class), 13
- bioAssayTuple,BioAssayTuples_assnlist-method (zzz-DTTCompatibility-class), 13
- BioAssayTuple-class (zzz-DTTCompatibility-class), 13
- bioAssayTuples_assnlist (zzz-DTTCompatibility-class), 13
- bioAssayTuples_assnlist,BioDataTuples-method (zzz-DTTCompatibility-class), 13
- BioAssayTuples_assnlist-class (zzz-DTTCompatibility-class), 13
- bioDataCube (zzz-DTTCompatibility-class), 13
- bioDataCube,BioDataValues_assn-method (zzz-DTTCompatibility-class), 13
- BioDataCube-class (zzz-DTTCompatibility-class), 13
- bioDataTuples (zzz-DTTCompatibility-class), 13
- bioDataTuples,BioDataValues_assn-method (zzz-DTTCompatibility-class), 13
- BioDataTuples-class (zzz-DTTCompatibility-class), 13
- bioDataValues_assn (zzz-DTTCompatibility-class), 13
- bioDataValues_assn,DerivedBioAssayData-method (zzz-DTTCompatibility-class), 13
- bioDataValues_assn,MeasuredBioAssayData-method (zzz-DTTCompatibility-class), 13
- BioDataValues_assn-class (zzz-DTTCompatibility-class), 13
- bioEvent_package (zzz-DTTCompatibility-class), 13
- bioEvent_package,MAGE-ML-method (zzz-DTTCompatibility-class), 13
- BioEvent_package-class (zzz-DTTCompatibility-class), 13
- biologicalCharacteristics_assnreflist (zzz-DTTCompatibility-class), 13
- biologicalCharacteristics_assnreflist,CompositeSequence (zzz-DTTCompatibility-class), 13
- BiologicalCharacteristics_assnreflist-class (zzz-DTTCompatibility-class), 13
- bioMaterial_assnlist (zzz-DTTCompatibility-class), 13
- bioMaterial_assnlist,BioMaterial_package-method (zzz-DTTCompatibility-class), 13

- 13
 BioMaterial_assnlist-class
 (zzz-DTTCompatibility-class),
 13
 bioMaterial_assnref
 (zzz-DTTCompatibility-class),
 13
 bioMaterial_assnref,BioMaterialMeasurement-method 13
 (zzz-DTTCompatibility-class),
 13
 bioMaterial_assnref,ManufactureLIMS-method
 (zzz-DTTCompatibility-class),
 13
 bioMaterial_assnref,ManufactureLIMSBiomaterial-method13
 (zzz-DTTCompatibility-class),
 13
 BioMaterial_assnref-class
 (zzz-DTTCompatibility-class),
 13
 bioMaterial_package
 (zzz-DTTCompatibility-class),
 13
 bioMaterial_package,MAGE-ML-method
 (zzz-DTTCompatibility-class),
 13
 BioMaterial_package-class
 (zzz-DTTCompatibility-class),
 13
 BioMaterial_ref-class
 (zzz-DTTCompatibility-class),
 13
 bioMaterialMeasurement
 (zzz-DTTCompatibility-class),
 13
 bioMaterialMeasurement,SourceBioMaterialMeasurements_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 BioMaterialMeasurement-class
 (zzz-DTTCompatibility-class),
 13
 bioMaterialPlateCol
 (zzz-DTTCompatibility-class),
 13
 bioMaterialPlateCol,ManufactureLIMSBiomaterial-method13
 (zzz-DTTCompatibility-class),
 13
 bioMaterialPlateIdentifier
 (zzz-DTTCompatibility-class),
 13
 bioMaterialPlateIdentifier,ManufactureLIMSBiomaterial-method
 (zzz-DTTCompatibility-class),
 13
 bioMaterialPlateRow
 (zzz-DTTCompatibility-class),
 13
 bioMaterialPlateRow,ManufactureLIMSBiomaterial-method
 (zzz-DTTCompatibility-class),
 13
 bioSample (zzz-DTTCompatibility-class),
 13
 bioSample,BioMaterial_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 BioSample-class
 (zzz-DTTCompatibility-class),
 13
 bioSample_ref
 (zzz-DTTCompatibility-class),
 13
 bioSample_ref,BioMaterial_assnref-method
 (zzz-DTTCompatibility-class),
 13
 BioSample_ref-class
 (zzz-DTTCompatibility-class),
 13
 bioSequence
 (zzz-DTTCompatibility-class),
 13
 bioSequence,BioSequence_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 BioSequence-class
 (zzz-DTTCompatibility-class),
 13
 bioSequence_assnlist
 (zzz-DTTCompatibility-class),
 13
 bioSequence_assnlist,BioSequence_package-method
 (zzz-DTTCompatibility-class),
 13
 BioSequence_assnlist-class
 (zzz-DTTCompatibility-class),
 13
 bioSequence_package
 (zzz-DTTCompatibility-class),
 13
 bioSequence_package,MAGE-ML-method
 (zzz-DTTCompatibility-class),
 13
 BioSequence_package-class
 (zzz-DTTCompatibility-class),
 13
 bioSequence_ref
 (zzz-DTTCompatibility-class),

- 13
- bioSequence_ref,BiologicalCharacteristics_assnreflist-method (zzz-DTTCompatibility-class), 13
- bioSequence_ref,ImmobilizedCharacteristics_assnreflist-method (zzz-DTTCompatibility-class), 13
- BioSequence_ref-class (zzz-DTTCompatibility-class), 13
- bioSource (zzz-DTTCompatibility-class), 13
- bioSource,BioMaterial_assnlist-method (zzz-DTTCompatibility-class), 13
- BioSource-class (zzz-DTTCompatibility-class), 13
- bioSource_ref (zzz-DTTCompatibility-class), 13
- bioSource_ref,BioMaterial_assnref-method (zzz-DTTCompatibility-class), 13
- BioSource_ref-class (zzz-DTTCompatibility-class), 13
- bqs_package (zzz-DTTCompatibility-class), 13
- bqs_package,MAGE-ML-method (zzz-DTTCompatibility-class), 13
- BQS_package-class (zzz-DTTCompatibility-class), 13
- category_assn (zzz-DTTCompatibility-class), 13
- category_assn,ExperimentalFactor-method (zzz-DTTCompatibility-class), 13
- Category_assn-class (zzz-DTTCompatibility-class), 13
- channel (zzz-DTTCompatibility-class), 13
- channel,Channel_assnlist-method (zzz-DTTCompatibility-class), 13
- Channel-class (zzz-DTTCompatibility-class), 13
- channel_assnlist (zzz-DTTCompatibility-class), 13
- channel_assnlist,BioAssay_package-method (zzz-DTTCompatibility-class), 13
- Channel_assnlist-class (zzz-DTTCompatibility-class), 13
- channel_assnref (zzz-DTTCompatibility-class), 13
- channel_assnref,DerivedSignal-method (zzz-DTTCompatibility-class), 13
- channel_assnref,Error-method (zzz-DTTCompatibility-class), 13
- channel_assnref,ExpectedValue-method (zzz-DTTCompatibility-class), 13
- channel_assnref,Failed-method (zzz-DTTCompatibility-class), 13
- channel_assnref,MeasuredSignal-method (zzz-DTTCompatibility-class), 13
- channel_assnref,PresentAbsent-method (zzz-DTTCompatibility-class), 13
- channel_assnref,PValue-method (zzz-DTTCompatibility-class), 13
- channel_assnref,Ratio-method (zzz-DTTCompatibility-class), 13
- channel_assnref,SpecializedQuantitationType-method (zzz-DTTCompatibility-class), 13
- Channel_assnref-class (zzz-DTTCompatibility-class), 13
- channel_ref (zzz-DTTCompatibility-class), 13
- channel_ref,Channel_assnref-method (zzz-DTTCompatibility-class), 13
- channel_ref,Channels_assnreflist-method (zzz-DTTCompatibility-class), 13
- Channel_ref-class

- (zzz-DTTCompatibility-class),
13
- channels_assnreflist
(zzz-DTTCompatibility-class),
13
- channels_assnreflist,DerivedBioAssay-method
(zzz-DTTCompatibility-class),
13
- channels_assnreflist,Image-method
(zzz-DTTCompatibility-class),
13
- channels_assnreflist,MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- channels_assnreflist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- Channels_assnreflist-class
(zzz-DTTCompatibility-class),
13
- characteristics_assnlist
(zzz-DTTCompatibility-class),
13
- characteristics_assnlist,BioSample-method
(zzz-DTTCompatibility-class),
13
- characteristics_assnlist,BioSource-method
(zzz-DTTCompatibility-class),
13
- characteristics_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- Characteristics_assnlist-class
(zzz-DTTCompatibility-class),
13
- clusterBioAssayData_assnref
(zzz-DTTCompatibility-class),
13
- clusterBioAssayData_assnref,BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
13
- ClusterBioAssayData_assnref-class
(zzz-DTTCompatibility-class),
13
- column (zzz-DTTCompatibility-class), 13
- column,FeatureLocation-method
(zzz-DTTCompatibility-class),
13
- column,Zone-method
(zzz-DTTCompatibility-class),
13
- componentCompounds_assnlist
(zzz-DTTCompatibility-class),
13
- componentCompounds_assnlist,Compound-method
(zzz-DTTCompatibility-class),
13
- ComponentCompounds_assnlist-class
(zzz-DTTCompatibility-class),
13
- composite_assnref
(zzz-DTTCompatibility-class),
13
- composite_assnref,CompositePosition-method
(zzz-DTTCompatibility-class),
13
- Composite_assnref-class
(zzz-DTTCompatibility-class),
13
- compositeCompositeMap
(zzz-DTTCompatibility-class),
13
- compositeCompositeMap,CompositeCompositeMap_assnlist-method
(zzz-DTTCompatibility-class),
13
- compositeCompositeMap,DesignElementMap_assnlist-method
(zzz-DTTCompatibility-class),
13
- CompositeCompositeMap-class
(zzz-DTTCompatibility-class),
13
- compositeCompositeMap_assnlist
(zzz-DTTCompatibility-class),
13
- compositeCompositeMap_assnlist,DesignElement_package-method
(zzz-DTTCompatibility-class),
13
- CompositeCompositeMap_assnlist-class
(zzz-DTTCompatibility-class),
13
- compositeCompositeMap_ref
(zzz-DTTCompatibility-class),
13
- compositeCompositeMap_ref,CompositeCompositeMaps_assnreflist
(zzz-DTTCompatibility-class),
13
- CompositeCompositeMap_ref-class
(zzz-DTTCompatibility-class),
13
- compositeCompositeMaps_assnreflist
(zzz-DTTCompatibility-class),

- 13
- compositeCompositeMaps_assnreflist, CompositeSequence-assnreflist-method (zzz-DTTCompatibility-class), 13
- CompositeCompositeMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- compositeGroup (zzz-DTTCompatibility-class), 13
- compositeGroup, CompositeGroup_assnlist-method (zzz-DTTCompatibility-class), 13
- CompositeGroup-class (zzz-DTTCompatibility-class), 13
- compositeGroup_assnlist (zzz-DTTCompatibility-class), 13
- compositeGroup_assnlist, ArrayDesign_package-method (zzz-DTTCompatibility-class), 13
- CompositeGroup_assnlist-class (zzz-DTTCompatibility-class), 13
- compositeGroup_ref (zzz-DTTCompatibility-class), 13
- compositeGroup_ref, CompositeGroups_assnreflist-method (zzz-DTTCompatibility-class), 13
- CompositeGroup_ref-class (zzz-DTTCompatibility-class), 13
- compositeGroups_assnreflist (zzz-DTTCompatibility-class), 13
- compositeGroups_assnreflist, ArrayDesign-method (zzz-DTTCompatibility-class), 13
- compositeGroups_assnreflist, PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 13
- CompositeGroups_assnreflist-class (zzz-DTTCompatibility-class), 13
- compositePosition (zzz-DTTCompatibility-class), 13
- compositePosition, CompositePositionSources_assnlist-method (zzz-DTTCompatibility-class), 13
- compositePosition, Coordinate_assn-method (zzz-DTTCompatibility-class), 13
- CompositePosition-class (zzz-DTTCompatibility-class), 13
- compositePositionSources_assnlist (zzz-DTTCompatibility-class), 13
- compositePositionSources_assnlist, CompositeCompositeMap-assnreflist-method (zzz-DTTCompatibility-class), 13
- CompositePositionSources_assnlist-class (zzz-DTTCompatibility-class), 13
- compositeSequence (zzz-DTTCompatibility-class), 13
- compositeSequence, CompositeSequence_assnlist-method (zzz-DTTCompatibility-class), 13
- compositeSequence, DesignElement_assnlist-method (zzz-DTTCompatibility-class), 13
- CompositeSequence-class (zzz-DTTCompatibility-class), 13
- compositeSequence_assnlist (zzz-DTTCompatibility-class), 13
- compositeSequence_assnlist, DesignElement_package-method (zzz-DTTCompatibility-class), 13
- CompositeSequence_assnlist-class (zzz-DTTCompatibility-class), 13
- compositeSequence_assnref (zzz-DTTCompatibility-class), 13
- compositeSequence_assnref, CompositeCompositeMap-method (zzz-DTTCompatibility-class), 13
- compositeSequence_assnref, ReporterCompositeMap-method (zzz-DTTCompatibility-class), 13
- CompositeSequence_assnref-class (zzz-DTTCompatibility-class), 13
- compositeSequence_ref (zzz-DTTCompatibility-class), 13
- compositeSequence_ref, Composite_assnref-method (zzz-DTTCompatibility-class), 13

- (zzz-DTTCompatibility-class),
13
- compound_assnlist
(zzz-DTTCompatibility-class),
13
- compound_assnlist, BioMaterial_package-method
(zzz-DTTCompatibility-class),
13
- compound_assnlist-class
(zzz-DTTCompatibility-class),
13
- compound_assnref
(zzz-DTTCompatibility-class),
13
- compound_assnref, CompoundMeasurement-method
(zzz-DTTCompatibility-class),
13
- Compound_assnref-class
(zzz-DTTCompatibility-class),
13
- compound_ref
(zzz-DTTCompatibility-class),
13
- compound_ref, Compound_assnref-method
(zzz-DTTCompatibility-class),
13
- compound_ref, Labels_assnreflist-method
(zzz-DTTCompatibility-class),
13
- Compound_ref-class
(zzz-DTTCompatibility-class),
13
- compoundIndices_assnlist
(zzz-DTTCompatibility-class),
13
- compoundIndices_assnlist, Compound-method
(zzz-DTTCompatibility-class),
13
- CompoundIndices_assnlist-class
(zzz-DTTCompatibility-class),
13
- compoundMeasurement
(zzz-DTTCompatibility-class),
13
- compoundMeasurement, ComponentCompounds_assnlist-method
(zzz-DTTCompatibility-class),
13
- compoundMeasurement, CompoundMeasurements_assnlist-method
(zzz-DTTCompatibility-class),
13
- CompoundMeasurement-class
(zzz-DTTCompatibility-class),
13
- compoundMeasurements_assnlist
- (zzz-DTTCompatibility-class),
13
- compositeSequence_ref, CompositeSequence_assnref-method
(zzz-DTTCompatibility-class),
13
- compositeSequence_ref, CompositeSequences_assnreflist-method
(zzz-DTTCompatibility-class),
13
- compositeSequence_ref, DesignElement_assnref-method
(zzz-DTTCompatibility-class),
13
- CompositeSequence_ref-class
(zzz-DTTCompatibility-class),
13
- compositeSequenceDimension
(zzz-DTTCompatibility-class),
13
- compositeSequenceDimension, DesignElementDimension_assnlist-method
(zzz-DTTCompatibility-class),
13
- CompositeSequenceDimension-class
(zzz-DTTCompatibility-class),
13
- compositeSequenceDimension_ref
(zzz-DTTCompatibility-class),
13
- compositeSequenceDimension_ref, DesignElementDimension_assnref-method
(zzz-DTTCompatibility-class),
13
- CompositeSequenceDimension_ref-class
(zzz-DTTCompatibility-class),
13
- compositeSequences_assnreflist
(zzz-DTTCompatibility-class),
13
- compositeSequences_assnreflist, CompositeGroup-method
(zzz-DTTCompatibility-class),
13
- compositeSequences_assnreflist, CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
13
- CompositeSequences_assnreflist-class
(zzz-DTTCompatibility-class),
13
- compound (zzz-DTTCompatibility-class),
13
- compound, Compound_assnlist-method
(zzz-DTTCompatibility-class),
13
- Compound-class
(zzz-DTTCompatibility-class),
13

- (zzz-DTTCompatibility-class),
13
- compoundMeasurements_assnlist, Treatment-method
(zzz-DTTCompatibility-class),
13
- CompoundMeasurements_assnlist-class
(zzz-DTTCompatibility-class),
13
- concentrationUnit
(zzz-DTTCompatibility-class),
13
- concentrationUnit, Unit_assn-method
(zzz-DTTCompatibility-class),
13
- ConcentrationUnit-class
(zzz-DTTCompatibility-class),
13
- ConfidenceIndicator_assnlist-class
(zzz-DTTCompatibility-class),
13
- ConfidenceIndicator_ref-class
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, DerivedSignal-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, Error-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, ExpectedValue-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, Failed-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, PValue-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, Ratio-method
(zzz-DTTCompatibility-class),
13
- confidenceIndicators_assnreflist, SpecializedQuantitiesType-assnreflist
(zzz-DTTCompatibility-class),
13
- ConfidenceIndicators_assnreflist-class
(zzz-DTTCompatibility-class),
13
- contact_assnlist
(zzz-DTTCompatibility-class),
13
- contact_assnlist, AuditAndSecurity_package-method
(zzz-DTTCompatibility-class),
13
- Contact_assnlist-class
(zzz-DTTCompatibility-class),
13
- Contact_ref-class
(zzz-DTTCompatibility-class),
13
- contacts_assnreflist
(zzz-DTTCompatibility-class),
13
- contacts_assnreflist, Database-method
(zzz-DTTCompatibility-class),
13
- Contacts_assnreflist-class
(zzz-DTTCompatibility-class),
13
- containedFeatures_assnreflist
(zzz-DTTCompatibility-class),
13
- containedFeatures_assnreflist, FeatureDimension-method
(zzz-DTTCompatibility-class),
13
- ContainedFeatures_assnreflist-class
(zzz-DTTCompatibility-class),
13
- control (zzz-ARRCompatibility-class), 12
- control, UserAttribute-method
(zzz-ARRCompatibility-class),
12
- Control-class
(zzz-ARRCompatibility-class),
12
- ControlFeatures_assnreflist
(zzz-DTTCompatibility-class),
13
- ControlFeatures_assnreflist, Feature-method
(zzz-DTTCompatibility-class),
13
- ControlFeatures_assnreflist-class
(zzz-DTTCompatibility-class),
13
- ControlFeatures_assnreflist, FeatureType-assnreflist
(zzz-DTTCompatibility-class),
13

- 13
- controlledFeatures_assnreflist, Feature-method
(zzz-DTTCompatibility-class),
13
- ControlledFeatures_assnreflist-class
(zzz-DTTCompatibility-class),
13
- controlType_assn
(zzz-DTTCompatibility-class),
13
- controlType_assn, CompositeSequence-method
(zzz-DTTCompatibility-class),
13
- controlType_assn, Feature-method
(zzz-DTTCompatibility-class),
13
- controlType_assn, Reporter-method
(zzz-DTTCompatibility-class),
13
- ControlType_assn-class
(zzz-DTTCompatibility-class),
13
- coordinate_assn
(zzz-DTTCompatibility-class),
13
- coordinate_assn, SeqFeatureLocation-method
(zzz-DTTCompatibility-class),
13
- Coordinate_assn-class
(zzz-DTTCompatibility-class),
13
- createdBy (zzz-ARRCompatibility-class),
12
- createdBy, ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- createdBy, PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- createdBy, TemplateFile-method
(zzz-ARRCompatibility-class),
12
- createdDateTime
(zzz-ARRCompatibility-class),
12
- createdDateTime, ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- createdDateTime, PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- createdDateTime, TemplateFile-method
(zzz-ARRCompatibility-class),
12
- createdStep
(zzz-ARRCompatibility-class),
12
- createdStep, ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- createdStep, PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- database (zzz-DTTCompatibility-class),
13
- database, Database_assnlist-method
(zzz-DTTCompatibility-class),
13
- Database-class
(zzz-DTTCompatibility-class),
13
- database_assnlist
(zzz-DTTCompatibility-class),
13
- database_assnlist, Description_package-method
(zzz-DTTCompatibility-class),
13
- Database_assnlist-class
(zzz-DTTCompatibility-class),
13
- database_assnref
(zzz-DTTCompatibility-class),
13
- database_assnref, DatabaseEntry-method
(zzz-DTTCompatibility-class),
13
- Database_assnref-class
(zzz-DTTCompatibility-class),
13
- database_ref
(zzz-DTTCompatibility-class),
13
- database_ref, Database_assnref-method
(zzz-DTTCompatibility-class),
13
- Database_ref-class
(zzz-DTTCompatibility-class),
13
- databaseEntry
(zzz-DTTCompatibility-class),
13
- databaseEntry, Accessions_assnlist-method
(zzz-DTTCompatibility-class),
13

- databaseEntry, DatabaseReferences_assnlist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class),
[13](#)
- databaseEntry, ExternalLIMS_assn-method
 (zzz-DTTCompatibility-class),
[13](#)
- databaseEntry, IdentifierLIMS_assn-method
 (zzz-DTTCompatibility-class),
[13](#)
- databaseEntry, OntologyReference_assn-method
 (zzz-DTTCompatibility-class),
[13](#)
- databaseEntry, SequenceDatabases_assnlist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class),
[13](#)
- DatabaseEntry-class
 (zzz-DTTCompatibility-class),
[13](#)
- databaseReferences_assnlist
 (zzz-DTTCompatibility-class),
[13](#)
- databaseReferences_assnlist, Description-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class),
[13](#)
- DatabaseReferences_assnlist-class
 (zzz-DTTCompatibility-class),
[13](#)
- dataExternal
 (zzz-DTTCompatibility-class),
[13](#)
- dataExternal, DataExternal_assn-method
 (zzz-DTTCompatibility-class),
[13](#)
- DataExternal-class
 (zzz-DTTCompatibility-class),
[13](#)
- dataExternal_assn
 (zzz-DTTCompatibility-class),
[13](#)
- dataExternal_assn, BioDataCube-method
 (zzz-DTTCompatibility-class),
[13](#)
- DataExternal_assn-class
 (zzz-DTTCompatibility-class),
[13](#)
- dataFormat
 (zzz-DTTCompatibility-class),
[13](#)
- dataFormat, DataExternal-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataFormatInfoURI
 (zzz-DTTCompatibility-class),
[13](#)
- dataFormatInfoURI, DataExternal-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataInternal
 (zzz-DTTCompatibility-class),
[13](#)
- dataInternal, DataInternal_assn-method
 (zzz-DTTCompatibility-class),
[13](#)
- DataInternal-class
 (zzz-DTTCompatibility-class),
[13](#)
- dataInternal_assn
 (zzz-DTTCompatibility-class),
[13](#)
- dataInternal_assn, BioDataCube-method
 (zzz-DTTCompatibility-class),
[13](#)
- DataInternal_assn-class
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, DerivedSignal-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, Error-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, ExpectedValue-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, Failed-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, MeasuredSignal-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, NodeValue-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, Parameter-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, PresentAbsent-method
 (zzz-DTTCompatibility-class),
[13](#)
- dataType_assn, PValue-method
 (zzz-DTTCompatibility-class),
[13](#)

- 13
- dataType_assn,Ratio-method
(zzz-DTTCompatibility-class),
13
- dataType_assn,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- DataType_assn-class
(zzz-DTTCompatibility-class),
13
- datum (zzz-DTTCompatibility-class), 13
- datum,Datum_assn-method
(zzz-DTTCompatibility-class),
13
- Datum-class
(zzz-DTTCompatibility-class),
13
- datum_assn
(zzz-DTTCompatibility-class),
13
- datum_assn,QuantitationTypeTuple-method
(zzz-DTTCompatibility-class),
13
- Datum_assn-class
(zzz-DTTCompatibility-class),
13
- defaultValue
(zzz-ARRCompatibility-class),
12
- defaultValue,UserAttribute-method
(zzz-ARRCompatibility-class),
12
- defaultValue_assn
(zzz-DTTCompatibility-class),
13
- defaultValue_assn,Parameter-method
(zzz-DTTCompatibility-class),
13
- DefaultValue_assn-class
(zzz-DTTCompatibility-class),
13
- defectType_assn
(zzz-DTTCompatibility-class),
13
- defectType_assn,FeatureDefect-method
(zzz-DTTCompatibility-class),
13
- defectType_assn,ZoneDefect-method
(zzz-DTTCompatibility-class),
13
- DefectType_assn-class
(zzz-DTTCompatibility-class),
13
- deltaX (zzz-DTTCompatibility-class), 13
- deltaX,PositionDelta-method
(zzz-DTTCompatibility-class),
13
- deltaY (zzz-DTTCompatibility-class), 13
- deltaY,PositionDelta-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssay
(zzz-DTTCompatibility-class),
13
- derivedBioAssay,BioAssay_assnlist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssay,DerivedBioAssay_assnlist-method
(zzz-DTTCompatibility-class),
13
- DerivedBioAssay-class
(zzz-DTTCompatibility-class),
13
- DerivedBioAssay_assnlist-class
(zzz-DTTCompatibility-class),
13
- derivedBioAssay_ref
(zzz-DTTCompatibility-class),
13
- derivedBioAssay_ref,BioAssay_assnref-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssay_ref,BioAssayMapTarget_assnref-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssay_ref,BioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssay_ref,SourceBioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssay_ref,TopLevelBioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
13
- DerivedBioAssay_ref-class
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData,BioAssayData_assnlist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData,DerivedBioAssayData_assnlist-method

- (zzz-DTTCompatibility-class),
13
- DerivedBioAssayData-class
(zzz-DTTCompatibility-class),
13
- DerivedBioAssayData_assnlist-class
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_assnreflist
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_assnreflist,DerivedBioAssay-method
(zzz-DTTCompatibility-class),
13
- DerivedBioAssayData_assnreflist-class
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_ref
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_ref,BioAssayData_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_ref,BioAssayDataSources_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_ref,ClusterBioAssayData_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_ref,DerivedBioAssayData_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedBioAssayData_ref,DerivedBioAssayDataTarget_assnreflist-method
(zzz-DTTCompatibility-class),
13
- DerivedBioAssayData_ref-class
(zzz-DTTCompatibility-class),
13
- DerivedBioAssayDataTarget_assnref-class
(zzz-DTTCompatibility-class),
13
- derivedBioAssayMap_assnreflist
(zzz-DTTCompatibility-class),
13
- derivedBioAssayMap_assnreflist,DerivedBioAssay-method
(zzz-DTTCompatibility-class),
13
- DerivedBioAssayMap_assnreflist-class
(zzz-DTTCompatibility-class),
13
- derivedSignal
(zzz-DTTCompatibility-class),
13
- derivedSignal,QuantitationType_assnlist-method
(zzz-DTTCompatibility-class),
13
- DerivedSignal-class
(zzz-DTTCompatibility-class),
13
- derivedSignal_ref
(zzz-DTTCompatibility-class),
13
- derivedSignal_ref,QuantitationType_assnref-method
(zzz-DTTCompatibility-class),
13
- derivedSignal_ref,QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedSignal_ref,SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class),
13
- derivedSignal_ref,TargetQuantitationType_assnref-method
(zzz-DTTCompatibility-class),
13
- DerivedSignal_ref-class
(zzz-DTTCompatibility-class),
13
- description
(zzz-DTTCompatibility-class),
13
- description,Descriptions_assnlist-method
(zzz-DTTCompatibility-class),
13
- description,NormalizationDescription_assn-method
(zzz-DTTCompatibility-class),
13
- description,OntologyEntry-method
(zzz-DTTCompatibility-class),
13
- description,QualityControlDescription_assn-method
(zzz-DTTCompatibility-class),
13
- description,ReplicateDescription_assn-method
(zzz-DTTCompatibility-class),
13
- Description-class
(zzz-DTTCompatibility-class),
13
- description_package
(zzz-DTTCompatibility-class),
13
- description_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
13

- Description_package-class (zzz-DTTCompatibility-class), 13
- descriptions_assnlist (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Array-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ArrayDesign-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ArrayGroup-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ArrayManufacture-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Audit-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BibliographicReference-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayCreation-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayDataCluster-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayDimension-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayMap-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayTreatment-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioSample-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioSequence-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioSource-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Channel-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeCompositeMap-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeGroup-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeSequence-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeSequenceDimension-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Compound-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Database-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,DerivedBioAssay-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,DerivedBioAssayData-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,DerivedSignal-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Description-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Error-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ExpectedValue-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Experiment-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ExperimentalFactor-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ExperimentDesign-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,FactorValue-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Failed-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Feature-method (zzz-DTTCompatibility-class), 13

- 13
- descriptions_assnlist, FeatureDimension-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, FeatureReporterMap-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Fiducial-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Hardware-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, HardwareApplication-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Hybridization-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Image-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, ManufactureLIMS-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Node-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, NodeContents-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Organization-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Parameter-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Person-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Protocol-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, PValue-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Ratio-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Reporter-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, ReporterDimension-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, ReporterGroup-method

- (zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Security-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, SecurityGroup-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, SeqFeature-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Software-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, SoftwareApplication-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Transformation-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Treatment-method
(zzz-DTTCompatibility-class),
13
- descriptions_assnlist, Zone-method
(zzz-DTTCompatibility-class),
13
- Descriptions_assnlist-class
(zzz-DTTCompatibility-class),
13
- DesignElement_assnlist-class
(zzz-DTTCompatibility-class),
13
- designElement_assnref
(zzz-DTTCompatibility-class),
13
- designElement_assnref, DesignElementTuple-method
(zzz-DTTCompatibility-class),
13
- DesignElement_assnref-class
(zzz-DTTCompatibility-class),
13
- designElement_package
(zzz-DTTCompatibility-class),
13
- designElement_package, MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- DesignElement_package-class
(zzz-DTTCompatibility-class),
13
- DesignElement_ref-class
(zzz-DTTCompatibility-class),
13
- designElementDimension_assnlist
(zzz-DTTCompatibility-class),
13
- designElementDimension_assnlist, BioAssayData_package-method
(zzz-DTTCompatibility-class),
13
- DesignElementDimension_assnlist-class
(zzz-DTTCompatibility-class),
13
- designElementDimension_assnref
(zzz-DTTCompatibility-class),
13
- designElementDimension_assnref, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
13
- designElementDimension_assnref, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- designElementDimension_assnref, NodeContents-method
(zzz-DTTCompatibility-class),
13
- DesignElementDimension_assnref-class
(zzz-DTTCompatibility-class),
13
- DesignElementDimension_ref-class
(zzz-DTTCompatibility-class),
13
- DesignElementMap_assnlist-class
(zzz-DTTCompatibility-class),
13
- DesignElementMap_ref-class
(zzz-DTTCompatibility-class),
13
- designElementMapping
(zzz-DTTCompatibility-class),
13
- designElementMapping, DesignElementMapping_assn-method
(zzz-DTTCompatibility-class),
13
- DesignElementMapping-class
(zzz-DTTCompatibility-class),
13
- designElementMapping_assn
(zzz-DTTCompatibility-class),
13
- designElementMapping_assn, Transformation-method
(zzz-DTTCompatibility-class),
13

- DesignElementMapping_assn-class (zzz-DTTCompatibility-class), 13
- DesignElementMapping_assn-class, (zzz-DTTCompatibility-class), 13
- designElementMaps_assnreflist (zzz-DTTCompatibility-class), 13
- designElementMaps_assnreflist, DesignElementMapping-method (zzz-DTTCompatibility-class), 13
- DesignElementMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- designElementTuple (zzz-DTTCompatibility-class), 13
- designElementTuple, DesignElementTuples_assnlist-method (zzz-DTTCompatibility-class), 13
- DesignElementTuple-class (zzz-DTTCompatibility-class), 13
- designElementTuples_assnlist (zzz-DTTCompatibility-class), 13
- designElementTuples_assnlist, BioAssayTuple-method (zzz-DTTCompatibility-class), 13
- DesignElementTuples_assnlist-class (zzz-DTTCompatibility-class), 13
- designProviders_assnreflist (zzz-DTTCompatibility-class), 13
- designProviders_assnreflist, ArrayDesign-method (zzz-DTTCompatibility-class), 13
- designProviders_assnreflist, PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 13
- DesignProviders_assnreflist-class (zzz-DTTCompatibility-class), 13
- directory (accessors), 2
- directory, NetAffxResource-method (NetAffxResource-class), 5
- distanceUnit (zzz-DTTCompatibility-class), 13
- distanceUnit, DistanceUnit_assn-method (zzz-DTTCompatibility-class), 13
- distanceUnit, Unit_assn-method (zzz-DTTCompatibility-class), 13
- DistanceUnit-class (zzz-DTTCompatibility-class), 13
- DistanceUnit-class (zzz-DTTCompatibility-class), 13
- distanceUnit_assn (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, ArrayGroup-method (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, FeatureGroup-method (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, Fiducial-method (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, Position-method (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, PositionDelta-method (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, Zone-method (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, ZoneGroup-method (zzz-DTTCompatibility-class), 13
- distanceUnit_assn, ZoneLayout-method (zzz-DTTCompatibility-class), 13
- DistanceUnit_assn-class (zzz-DTTCompatibility-class), 13
- DTTCompatibility, 13
- DTTCompatibility-class (AffyCompatible-class), 3
- editor (zzz-DTTCompatibility-class), 13
- editor, BibliographicReference-method (zzz-DTTCompatibility-class), 13
- email (zzz-DTTCompatibility-class), 13
- email, Organization-method (zzz-DTTCompatibility-class), 13
- email, Person-method (zzz-DTTCompatibility-class), 13
- error (zzz-DTTCompatibility-class), 13
- error, ConfidenceIndicator_assnlist-method (zzz-DTTCompatibility-class), 13

13
 error, QuantitationType_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 Error-class
 (zzz-DTTCompatibility-class),
 13
 error_ref (zzz-DTTCompatibility-class),
 13
 error_ref, ConfidenceIndicators_assnreflist-method
 (zzz-DTTCompatibility-class),
 13
 error_ref, QuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 13
 error_ref, QuantitationTypes_assnreflist-method
 (zzz-DTTCompatibility-class),
 13
 error_ref, SourcesQuantitationType_assnreflist-method
 (zzz-DTTCompatibility-class),
 13
 error_ref, TargetQuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 13
 Error_ref-class
 (zzz-DTTCompatibility-class),
 13
 expectedValue
 (zzz-DTTCompatibility-class),
 13
 expectedValue, ConfidenceIndicator_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 expectedValue, QuantitationType_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 ExpectedValue-class
 (zzz-DTTCompatibility-class),
 13
 expectedValue_ref
 (zzz-DTTCompatibility-class),
 13
 expectedValue_ref, ConfidenceIndicators_assnreflist-method
 (zzz-DTTCompatibility-class),
 13
 expectedValue_ref, QuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 13
 expectedValue_ref, QuantitationTypes_assnreflist-method
 (zzz-DTTCompatibility-class),
 13
 expectedValue_ref, SourcesQuantitationType_assnreflist-method
 (zzz-DTTCompatibility-class),
 13
 (zzz-DTTCompatibility-class),
 13
 expectedValue_ref, TargetQuantitationType_assnref-method
 (zzz-DTTCompatibility-class),
 13
 ExpectedValue_ref-class
 (zzz-DTTCompatibility-class),
 13
 experiment
 (zzz-DTTCompatibility-class),
 13
 experiment, Experiment_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 Experiment-class
 (zzz-DTTCompatibility-class),
 13
 experiment_assnlist
 (zzz-DTTCompatibility-class),
 13
 experiment_assnlist, Experiment_package-method
 (zzz-DTTCompatibility-class),
 13
 Experiment_assnlist-class
 (zzz-DTTCompatibility-class),
 13
 experiment_package
 (zzz-DTTCompatibility-class),
 13
 experiment_package, MAGE-ML-method
 (zzz-DTTCompatibility-class),
 13
 Experiment_package-class
 (zzz-DTTCompatibility-class),
 13
 experimentalFactor
 (zzz-DTTCompatibility-class),
 13
 experimentalFactor, ExperimentalFactors_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 ExperimentalFactor-class
 (zzz-DTTCompatibility-class),
 13
 ExperimentalFactor_assnref-class
 (zzz-DTTCompatibility-class),
 13
 experimentalFactor_ref
 (zzz-DTTCompatibility-class),
 13
 experimentalFactor_ref, ExperimentalFactor_assnref-method
 (zzz-DTTCompatibility-class),
 13

- 13
 ExperimentalFactor_ref-class
 (zzz-DTTCompatibility-class),
 13
 experimentalFactors_assnlist
 (zzz-DTTCompatibility-class),
 13
 experimentalFactors_assnlist,ExperimentDesign-method 13
 (zzz-DTTCompatibility-class),
 13
 ExperimentalFactors_assnlist-class
 (zzz-DTTCompatibility-class),
 13
 experimentDesign
 (zzz-DTTCompatibility-class),
 13
 experimentDesign,ExperimentDesigns_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 ExperimentDesign-class
 (zzz-DTTCompatibility-class),
 13
 experimentDesigns_assnlist
 (zzz-DTTCompatibility-class),
 13
 experimentDesigns_assnlist,Experiment-method
 (zzz-DTTCompatibility-class),
 13
 ExperimentDesigns_assnlist-class
 (zzz-DTTCompatibility-class),
 13
 exportedFromDB
 (zzz-DTTCompatibility-class),
 13
 exportedFromDB,ExternalReference-method
 (zzz-DTTCompatibility-class),
 13
 exportedFromServer
 (zzz-DTTCompatibility-class),
 13
 exportedFromServer,ExternalReference-method
 (zzz-DTTCompatibility-class),
 13
 exportID (zzz-DTTCompatibility-class),
 13
 exportID,ExternalReference-method
 (zzz-DTTCompatibility-class),
 13
 exportName
 (zzz-DTTCompatibility-class),
 13
 exportName,ExternalReference-method
 (zzz-DTTCompatibility-class),
 13
 externalLIMS_assn
 (zzz-DTTCompatibility-class),
 13
 externalLIMS_assn,Compound-method
 (zzz-DTTCompatibility-class),
 13
 ExternalLIMS_assn-class
 (zzz-DTTCompatibility-class),
 13
 externalReference
 (zzz-DTTCompatibility-class),
 13
 externalReference,ExternalReference_assn-method
 (zzz-DTTCompatibility-class),
 13
 externalReference_assn
 (zzz-DTTCompatibility-class),
 13
 externalReference_assn,Description-method
 (zzz-DTTCompatibility-class),
 13
 ExternalReference_assn-class
 (zzz-DTTCompatibility-class),
 13
 factorValue
 (zzz-DTTCompatibility-class),
 13
 factorValue,FactorValues_assnlist-method
 (zzz-DTTCompatibility-class),
 13
 FactorValue-class
 (zzz-DTTCompatibility-class),
 13
 factorValue_ref
 (zzz-DTTCompatibility-class),
 13
 factorValue_ref,BioAssayFactorValues_assnreflist-method
 (zzz-DTTCompatibility-class),
 13
 FactorValue_ref-class
 (zzz-DTTCompatibility-class),
 13
 factorValues_assnlist
 (zzz-DTTCompatibility-class),
 13
 factorValues_assnlist,ExperimentalFactor-method
 (zzz-DTTCompatibility-class),
 13
 FactorValues_assnlist-class
 (zzz-DTTCompatibility-class),

- 13
- failed (zzz-DTTCompatibility-class), 13
- failed, QuantitationType_assnlist-method (zzz-DTTCompatibility-class), 13
- Failed-class (zzz-DTTCompatibility-class), 13
- failed_ref (zzz-DTTCompatibility-class), 13
- failed_ref, QuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- failed_ref, QuantitationTypes_assnreflist-method (zzz-DTTCompatibility-class), 13
- failed_ref, SourcesQuantitationType_assnreflist-method (zzz-DTTCompatibility-class), 13
- failed_ref, TargetQuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- Failed_ref-class (zzz-DTTCompatibility-class), 13
- failTypes_assnlist (zzz-DTTCompatibility-class), 13
- failTypes_assnlist, Reporter-method (zzz-DTTCompatibility-class), 13
- FailTypes_assnlist-class (zzz-DTTCompatibility-class), 13
- fax (zzz-DTTCompatibility-class), 13
- fax, Organization-method (zzz-DTTCompatibility-class), 13
- fax, Person-method (zzz-DTTCompatibility-class), 13
- feature (zzz-DTTCompatibility-class), 13
- feature, DesignElement_assnlist-method (zzz-DTTCompatibility-class), 13
- feature, Features_assnlist-method (zzz-DTTCompatibility-class), 13
- Feature-class (zzz-DTTCompatibility-class), 13
- feature_assnref (zzz-DTTCompatibility-class), 13
- feature_assnref, FeatureDefect-method (zzz-DTTCompatibility-class), 13
- feature_assnref, FeatureInformation-method (zzz-DTTCompatibility-class), 13
- feature_assnref, ManufactureLIMS-method (zzz-DTTCompatibility-class), 13
- feature_assnref, ManufactureLIMSBiomaterial-method (zzz-DTTCompatibility-class), 13
- Feature_assnref-class (zzz-DTTCompatibility-class), 13
- feature_ref (zzz-DTTCompatibility-class), 13
- feature_ref, ContainedFeatures_assnreflist-method (zzz-DTTCompatibility-class), 13
- feature_ref, ControlFeatures_assnreflist-method (zzz-DTTCompatibility-class), 13
- feature_ref, ControlledFeatures_assnreflist-method (zzz-DTTCompatibility-class), 13
- feature_ref, DesignElement_assnref-method (zzz-DTTCompatibility-class), 13
- feature_ref, Feature_assnref-method (zzz-DTTCompatibility-class), 13
- Feature_ref-class (zzz-DTTCompatibility-class), 13
- featureDefect (zzz-DTTCompatibility-class), 13
- featureDefect, FeatureDefects_assnlist-method (zzz-DTTCompatibility-class), 13
- FeatureDefect-class (zzz-DTTCompatibility-class), 13
- featureDefects_assnlist (zzz-DTTCompatibility-class), 13
- featureDefects_assnlist, ArrayManufactureDeviation-method

- (zzz-DTTCompatibility-class),
13
- FeatureDefects_assnlist-class
(zzz-DTTCompatibility-class),
13
- featureDimension
(zzz-DTTCompatibility-class),
13
- featureDimension,DesignElementDimension_assnlist-method
(zzz-DTTCompatibility-class),
13
- FeatureDimension-class
(zzz-DTTCompatibility-class),
13
- featureDimension_ref
(zzz-DTTCompatibility-class),
13
- featureDimension_ref,DesignElementDimension_assnref-method
(zzz-DTTCompatibility-class),
13
- FeatureDimension_ref-class
(zzz-DTTCompatibility-class),
13
- featureExtraction
(zzz-DTTCompatibility-class),
13
- featureExtraction,FeatureExtraction_assn-method
(zzz-DTTCompatibility-class),
13
- FeatureExtraction-class
(zzz-DTTCompatibility-class),
13
- featureExtraction_assn
(zzz-DTTCompatibility-class),
13
- featureExtraction_assn,MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- FeatureExtraction_assn-class
(zzz-DTTCompatibility-class),
13
- featureGroup
(zzz-DTTCompatibility-class),
13
- featureGroup,FeatureGroups_assnlist-method
(zzz-DTTCompatibility-class),
13
- FeatureGroup-class
(zzz-DTTCompatibility-class),
13
- FeatureGroup_assnref-class
(zzz-DTTCompatibility-class),
13
- featureGroup_ref
(zzz-DTTCompatibility-class),
13
- featureGroup_ref,FeatureGroup_assnref-method
(zzz-DTTCompatibility-class),
13
- FeatureGroup_ref-class
(zzz-DTTCompatibility-class),
13
- featureGroups_assnlist
(zzz-DTTCompatibility-class),
13
- featureGroups_assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class),
13
- featureGroups_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- FeatureGroups_assnlist-class
(zzz-DTTCompatibility-class),
13
- featureHeight
(zzz-DTTCompatibility-class),
13
- featureHeight,FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- featureInformation
(zzz-DTTCompatibility-class),
13
- featureInformation,FeatureInformationSources_assnlist-method
(zzz-DTTCompatibility-class),
13
- FeatureInformation-class
(zzz-DTTCompatibility-class),
13
- featureInformationSources_assnlist
(zzz-DTTCompatibility-class),
13
- featureInformationSources_assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class),
13
- FeatureInformationSources_assnlist-class
(zzz-DTTCompatibility-class),
13
- featureLength
(zzz-DTTCompatibility-class),
13
- featureLength,FeatureGroup-method
(zzz-DTTCompatibility-class),
13

- featureIMSSs_assnlist (zzz-DTTCompatibility-class), 13
- featureIMSSs_assnlist,ArrayManufacture-method (zzz-DTTCompatibility-class), 13
- FeatureIMSSs_assnlist-class (zzz-DTTCompatibility-class), 13
- featureLocation (zzz-DTTCompatibility-class), 13
- featureLocation,FeatureLocation_assn-method (zzz-DTTCompatibility-class), 13
- FeatureLocation-class (zzz-DTTCompatibility-class), 13
- featureLocation_assn (zzz-DTTCompatibility-class), 13
- featureLocation_assn,Feature-method (zzz-DTTCompatibility-class), 13
- FeatureLocation_assn-class (zzz-DTTCompatibility-class), 13
- featureReporterMap (zzz-DTTCompatibility-class), 13
- featureReporterMap,DesignElementMap_assnlist-method (zzz-DTTCompatibility-class), 13
- featureReporterMap,FeatureReporterMap_assnlist-method (zzz-DTTCompatibility-class), 13
- FeatureReporterMap-class (zzz-DTTCompatibility-class), 13
- featureReporterMap_assnlist (zzz-DTTCompatibility-class), 13
- featureReporterMap_assnlist,DesignElement_package-method (zzz-DTTCompatibility-class), 13
- FeatureReporterMap_assnlist-class (zzz-DTTCompatibility-class), 13
- featureReporterMap_ref (zzz-DTTCompatibility-class), 13
- featureReporterMap_ref,DesignElementMaps_assnreflist-method (zzz-DTTCompatibility-class), 13
- featureReporterMap_ref,FeatureReporterMaps_assnreflist-method (zzz-DTTCompatibility-class), 13
- FeatureReporterMap_ref-class (zzz-DTTCompatibility-class), 13
- featureReporterMaps_assnreflist (zzz-DTTCompatibility-class), 13
- featureReporterMaps_assnreflist,Reporter-method (zzz-DTTCompatibility-class), 13
- FeatureReporterMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- features_assnlist (zzz-DTTCompatibility-class), 13
- features_assnlist,FeatureGroup-method (zzz-DTTCompatibility-class), 13
- Features_assnlist-class (zzz-DTTCompatibility-class), 13
- featureShape_assn (zzz-DTTCompatibility-class), 13
- featureShape_assn,FeatureGroup-method (zzz-DTTCompatibility-class), 13
- FeatureShape_assn-class (zzz-DTTCompatibility-class), 13
- featureWidth (zzz-DTTCompatibility-class), 13
- featureWidth,FeatureGroup-method (zzz-DTTCompatibility-class), 13
- fiducial (zzz-DTTCompatibility-class), 13
- fiducial,Fiducials_assnlist-method (zzz-DTTCompatibility-class), 13
- Fiducial-class (zzz-DTTCompatibility-class), 13
- fiducials_assnlist (zzz-DTTCompatibility-class), 13

- fiducials_assnlist, ArrayGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- Fiducials_assnlist-class
(zzz-DTTCompatibility-class),
[13](#)
- fiducialType_assn
(zzz-DTTCompatibility-class),
[13](#)
- fiducialType_assn, Fiducial-method
(zzz-DTTCompatibility-class),
[13](#)
- FiducialType_assn-class
(zzz-DTTCompatibility-class),
[13](#)
- filenameURI
(zzz-DTTCompatibility-class),
[13](#)
- filenameURI, DataExternal-method
(zzz-DTTCompatibility-class),
[13](#)
- firstName (zzz-DTTCompatibility-class),
[13](#)
- firstName, Person-method
(zzz-DTTCompatibility-class),
[13](#)
- format_assn
(zzz-DTTCompatibility-class),
[13](#)
- format_assn, Image-method
(zzz-DTTCompatibility-class),
[13](#)
- Format_assn-class
(zzz-DTTCompatibility-class),
[13](#)
- guid (zzz-ARRCompatibility-class), [12](#)
- guid, ArraySetFile-method
(zzz-ARRCompatibility-class),
[12](#)
- guid, PhysicalArray-method
(zzz-ARRCompatibility-class),
[12](#)
- guid, TemplateFile-method
(zzz-ARRCompatibility-class),
[12](#)
- hardware (zzz-DTTCompatibility-class),
[13](#)
- hardware, Hardware_assnlist-method
(zzz-DTTCompatibility-class),
[13](#)
- Hardware-class
(zzz-DTTCompatibility-class),
[13](#)
- hardware_assnlist
(zzz-DTTCompatibility-class),
[13](#)
- hardware_assnlist, Protocol_package-method
(zzz-DTTCompatibility-class),
[13](#)
- Hardware_assnlist-class
(zzz-DTTCompatibility-class),
[13](#)
- hardware_assnref
(zzz-DTTCompatibility-class),
[13](#)
- hardware_assnref, HardwareApplication-method
(zzz-DTTCompatibility-class),
[13](#)
- hardware_assnref, Software-method
(zzz-DTTCompatibility-class),
[13](#)
- Hardware_assnref-class
(zzz-DTTCompatibility-class),
[13](#)
- hardware_ref
(zzz-DTTCompatibility-class),
[13](#)
- hardware_ref, Hardware_assnref-method
(zzz-DTTCompatibility-class),
[13](#)
- hardware_ref, Hardwares_assnreflist-method
(zzz-DTTCompatibility-class),
[13](#)
- Hardware_ref-class
(zzz-DTTCompatibility-class),
[13](#)
- hardwareApplication
(zzz-DTTCompatibility-class),
[13](#)
- hardwareApplication, HardwareApplications_assnlist-method
(zzz-DTTCompatibility-class),
[13](#)
- HardwareApplication-class
(zzz-DTTCompatibility-class),
[13](#)
- hardwareApplications_assnlist
(zzz-DTTCompatibility-class),
[13](#)
- hardwareApplications_assnlist, ProtocolApplication-method
(zzz-DTTCompatibility-class),
[13](#)
- HardwareApplications_assnlist-class

- (zzz-DTTCompatibility-class),
13
- hardwareManufacturers_assnreflist
(zzz-DTTCompatibility-class),
13
- hardwareManufacturers_assnreflist,Hardware-method
(zzz-DTTCompatibility-class),
13
- HardwareManufacturers_assnreflist-class
(zzz-DTTCompatibility-class),
13
- hardwares_assnreflist
(zzz-DTTCompatibility-class),
13
- hardwares_assnreflist,Protocol-method
(zzz-DTTCompatibility-class),
13
- Hardwares_assnreflist-class
(zzz-DTTCompatibility-class),
13
- higherLevelAnalysis_package
(zzz-DTTCompatibility-class),
13
- higherLevelAnalysis_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- HigherLevelAnalysis_package-class
(zzz-DTTCompatibility-class),
13
- hybridization
(zzz-DTTCompatibility-class),
13
- hybridization,BioAssayCreation_assn-method
(zzz-DTTCompatibility-class),
13
- Hybridization-class
(zzz-DTTCompatibility-class),
13
- identifier
(zzz-DTTCompatibility-class),
13
- identifier,Array-method
(zzz-DTTCompatibility-class),
13
- identifier,Array_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,ArrayDesign-method
(zzz-DTTCompatibility-class),
13
- identifier,ArrayDesign_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- identifier,ArrayGroup_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,ArrayManufacture-method
(zzz-DTTCompatibility-class),
13
- identifier,ArrayManufacture_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssay_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayCreation-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayData_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayDataCluster_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayDimension-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayDimension_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayMap-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayMap_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,BioAssayTreatment-method
(zzz-DTTCompatibility-class),
13
- identifier,BioMaterial_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,BioSample-method
(zzz-DTTCompatibility-class),
13
- identifier,BioSample_ref-method
(zzz-DTTCompatibility-class),
13

- identifier, BioSequence-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, BioSequence_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, BioSource-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, BioSource_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Channel-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Channel_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeCompositeMap_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeGroup_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeSequence-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeSequence_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, CompositeSequenceDimension_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Compound-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Compound_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, ConfidenceIndicator_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Contact_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Database-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Database_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DerivedBioAssay_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DerivedBioAssayData_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DerivedSignal-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DerivedSignal_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DesignElement_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DesignElementDimension_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, DesignElementMap_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Error-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Error_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, ExpectedValue-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, ExpectedValue_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, Experiment-method
(zzz-DTTCompatibility-class),
[13](#)
- identifier, ExperimentalFactor-method
(zzz-DTTCompatibility-class),

- 13
- identifier, ExperimentalFactor_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, FactorValue-method
(zzz-DTTCompatibility-class),
13
- identifier, FactorValue_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, Failed-method
(zzz-DTTCompatibility-class),
13
- identifier, Failed_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, Feature-method
(zzz-DTTCompatibility-class),
13
- identifier, Feature_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, FeatureDimension-method
(zzz-DTTCompatibility-class),
13
- identifier, FeatureDimension_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- identifier, FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- identifier, FeatureGroup_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, FeatureReporterMap-method
(zzz-DTTCompatibility-class),
13
- identifier, FeatureReporterMap_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, Hardware-method
(zzz-DTTCompatibility-class),
13
- identifier, Hardware_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, Hybridization-method
(zzz-DTTCompatibility-class),
13
- identifier, Image-method
(zzz-DTTCompatibility-class),
13
- identifier, Image_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
- identifier, LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- identifier, LabeledExtract_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- identifier, MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- identifier, MeasuredBioAssay_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- identifier, MeasuredBioAssayData_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- identifier, MeasuredSignal_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, Organization-method
(zzz-DTTCompatibility-class),
13
- identifier, Organization_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, Parameter-method
(zzz-DTTCompatibility-class),
13
- identifier, Parameter_ref-method
(zzz-DTTCompatibility-class),
13
- identifier, Person-method
(zzz-DTTCompatibility-class),
13
- identifier, Person_ref-method

- (zzz-DTTCompatibility-class),
13
- identifier,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- identifier,PhysicalArrayDesign_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- identifier,PhysicalBioAssay_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- identifier,PresentAbsent_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,Protocol-method
(zzz-DTTCompatibility-class),
13
- identifier,Protocol_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,PValue-method
(zzz-DTTCompatibility-class),
13
- identifier,PValue_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,QuantitationType_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
13
- identifier,QuantitationTypeDimension_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
13
- identifier,QuantitationTypeMap_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,Ratio-method
(zzz-DTTCompatibility-class),
13
- identifier,Ratio_ref-method
(zzz-DTTCompatibility-class),
13
- 13
- identifier,Reporter-method
(zzz-DTTCompatibility-class),
13
- identifier,Reporter_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
13
- identifier,ReporterCompositeMap_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,ReporterDimension-method
(zzz-DTTCompatibility-class),
13
- identifier,ReporterDimension_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,ReporterGroup-method
(zzz-DTTCompatibility-class),
13
- identifier,ReporterGroup_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,Security-method
(zzz-DTTCompatibility-class),
13
- identifier,Security_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,SecurityGroup-method
(zzz-DTTCompatibility-class),
13
- identifier,SecurityGroup_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,Software-method
(zzz-DTTCompatibility-class),
13
- identifier,Software_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- identifier,SpecializedQuantitationType_ref-method
(zzz-DTTCompatibility-class),
13
- identifier,StandardQuantitationType_ref-method
(zzz-DTTCompatibility-class),
13

- identifieur,Transformation-method
(zzz-DTTCompatibility-class),
13
 - identifieur,Treatment-method
(zzz-DTTCompatibility-class),
13
 - identifieur,Zone-method
(zzz-DTTCompatibility-class),
13
 - identifieur,Zone_ref-method
(zzz-DTTCompatibility-class),
13
 - identifieurLIMS_assn
(zzz-DTTCompatibility-class),
13
 - identifieurLIMS_assn,ManufactureLIMS-method
(zzz-DTTCompatibility-class),
13
 - identifieurLIMS_assn,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
13
 - IdentifieurLIMS_assn-class
(zzz-DTTCompatibility-class),
13
 - Image-class
(zzz-DTTCompatibility-class),
13
 - image_ref (zzz-DTTCompatibility-class),
13
 - image_ref,Images_assnreflist-method
(zzz-DTTCompatibility-class),
13
 - Image_ref-class
(zzz-DTTCompatibility-class),
13
 - imageAcquisition
(zzz-DTTCompatibility-class),
13
 - imageAcquisition,BioAssayTreatments_assnreflist-method
(zzz-DTTCompatibility-class),
13
 - ImageAcquisition-class
(zzz-DTTCompatibility-class),
13
 - images_assnreflist
(zzz-DTTCompatibility-class),
13
 - images_assnreflist,ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
 - Images_assnreflist-class
(zzz-DTTCompatibility-class),
- 13
- immobilizedCharacteristics_assnreflist
(zzz-DTTCompatibility-class),
13
- immobilizedCharacteristics_assnreflist,Reporter-method
(zzz-DTTCompatibility-class),
13
- ImmobilizedCharacteristics_assnreflist-class
(zzz-DTTCompatibility-class),
13
- information_assnref
(zzz-DTTCompatibility-class),
13
- information_assnref,Array-method
(zzz-DTTCompatibility-class),
13
- Information_assnref-class
(zzz-DTTCompatibility-class),
- Internal-class, 4
- isApproximateLength
(zzz-DTTCompatibility-class),
13
- isApproximateLength,BioSequence-method
(zzz-DTTCompatibility-class),
13
- isBackground
(zzz-DTTCompatibility-class),
13
- isBackground,DerivedSignal-method
(zzz-DTTCompatibility-class),
13
- isBackground,Error-method
(zzz-DTTCompatibility-class),
13
- isBackground,ExpectedValue-method
(zzz-DTTCompatibility-class),
13
- isBackground,Failed-method
(zzz-DTTCompatibility-class),
13
- isBackground,MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- isBackground,PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- isBackground,PValue-method
(zzz-DTTCompatibility-class),
13
- isBackground,Ratio-method
(zzz-DTTCompatibility-class),

- 13
- isBackground, SpecializedQuantitationType-method (zzz-DTTCompatibility-class), 13
- isCircular (zzz-DTTCompatibility-class), 13
- isCircular, BioSequence-method (zzz-DTTCompatibility-class), 13
- isSolvent (zzz-DTTCompatibility-class), 13
- isSolvent, Compound-method (zzz-DTTCompatibility-class), 13
- issue (zzz-DTTCompatibility-class), 13
- issue, BibliographicReference-method (zzz-DTTCompatibility-class), 13
- kindCV (zzz-DTTCompatibility-class), 13
- kindCV, Measurement-method (zzz-DTTCompatibility-class), 13
- labeledExtract (zzz-DTTCompatibility-class), 13
- labeledExtract, BioMaterial_assnlist-method (zzz-DTTCompatibility-class), 13
- LabeledExtract-class (zzz-DTTCompatibility-class), 13
- labeledExtract_ref (zzz-DTTCompatibility-class), 13
- labeledExtract_ref, BioMaterial_assnref-method (zzz-DTTCompatibility-class), 13
- LabeledExtract_ref-class (zzz-DTTCompatibility-class), 13
- labels_assnreflist (zzz-DTTCompatibility-class), 13
- labels_assnreflist, Channel-method (zzz-DTTCompatibility-class), 13
- labels_assnreflist, LabeledExtract-method (zzz-DTTCompatibility-class), 13
- Labels_assnreflist-class (zzz-DTTCompatibility-class), 13
- lastName (zzz-DTTCompatibility-class), 13
- lastName, Person-method (zzz-DTTCompatibility-class), 13
- length, ArrayGroup-method (zzz-DTTCompatibility-class), 13
- length, BioSequence-method (zzz-DTTCompatibility-class), 13
- libraryPackageName (zzz-ARRCompatibility-class), 12
- libraryPackageName, PhysicalArray-method (zzz-ARRCompatibility-class), 12
- list, 4
- lowerRightX (zzz-DTTCompatibility-class), 13
- lowerRightX, Zone-method (zzz-DTTCompatibility-class), 13
- lowerRightY (zzz-DTTCompatibility-class), 13
- lowerRightY, Zone-method (zzz-DTTCompatibility-class), 13
- MAGE-ML-class (zzz-DTTCompatibility-class), 13
- make (zzz-DTTCompatibility-class), 13
- make, Hardware-method (zzz-DTTCompatibility-class), 13
- manufactureLIMS (zzz-DTTCompatibility-class), 13
- manufactureLIMS, FeatureLIMSS_assnlist-method (zzz-DTTCompatibility-class), 13
- ManufactureLIMS-class (zzz-DTTCompatibility-class), 13
- manufactureLIMSBiomaterial (zzz-DTTCompatibility-class), 13

- manufactureLIMSBiomaterial, FeatureLIMSS_assnlist-method
 (zzz-DTTCompatibility-class),
[13](#)
- ManufactureLIMSBiomaterial-class
 (zzz-DTTCompatibility-class),
[13](#)
- manufacturingDate
 (zzz-DTTCompatibility-class),
[13](#)
- manufacturingDate, ArrayManufacture-method
 (zzz-DTTCompatibility-class),
[13](#)
- massUnit (zzz-DTTCompatibility-class),
[13](#)
- massUnit, Unit_assn-method
 (zzz-DTTCompatibility-class),
[13](#)
- MassUnit-class
 (zzz-DTTCompatibility-class),
[13](#)
- masterFileGUID
 (zzz-ARRCompatibility-class),
[12](#)
- masterFileGUID, PhysicalArray-method
 (zzz-ARRCompatibility-class),
[12](#)
- masterFileName
 (zzz-ARRCompatibility-class),
[12](#)
- masterFileName, PhysicalArray-method
 (zzz-ARRCompatibility-class),
[12](#)
- materialType_assn
 (zzz-DTTCompatibility-class),
[13](#)
- materialType_assn, BioSample-method
 (zzz-DTTCompatibility-class),
[13](#)
- materialType_assn, BioSource-method
 (zzz-DTTCompatibility-class),
[13](#)
- materialType_assn, LabeledExtract-method
 (zzz-DTTCompatibility-class),
[13](#)
- MaterialType_assn-class
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay, BioAssay_assnlist-method
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay, MeasuredBioAssay_assnlist-method
 (zzz-DTTCompatibility-class),
[13](#)
- MeasuredBioAssay-class
 (zzz-DTTCompatibility-class),
[13](#)
- MeasuredBioAssay_assnlist-class
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay_ref
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay_ref, BioAssay_assnref-method
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay_ref, BioAssays_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay_ref, MeasuredBioAssayTarget_assnref-method
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay_ref, SourceBioAssays_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssay_ref, TopLevelBioAssays_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
- MeasuredBioAssay_ref-class
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssayData
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssayData, BioAssayData_assnlist-method
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssayData, MeasuredBioAssayData_assnlist-method
 (zzz-DTTCompatibility-class),
[13](#)
- MeasuredBioAssayData-class
 (zzz-DTTCompatibility-class),
[13](#)
- MeasuredBioAssayData_assnlist-class
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssayData_assnreflist
 (zzz-DTTCompatibility-class),
[13](#)
- measuredBioAssayData_assnreflist, MeasuredBioAssay-method
 (zzz-DTTCompatibility-class),
[13](#)

MeasuredBioAssayData_assnreflist-class (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement, ActionMeasurement_assn-method
 measuredBioAssayData_ref (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement, DefaultValue_assn-method
 measuredBioAssayData_ref, BioAssayData_assnreflist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement, Measurement_assn-method
 measuredBioAssayData_ref, BioAssayDataSources_assnreflist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 Measurement-class
 measuredBioAssayData_ref, ClusterBioAssayData_assnreflist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement_assn
 measuredBioAssayData_ref, MeasuredBioAssayData_assnreflist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement_assn, BioMaterialMeasurement-method
 MeasuredBioAssayData_ref-class (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement_assn, CompoundMeasurement-method
 MeasuredBioAssayTarget_assnref-class (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement_assn, FactorValue-method
 measuredSignal (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 Measurement_assn-class
 measuredSignal, QuantitationType_assnlist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement_package
 MeasuredSignal-class (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 measurement_package, MAGE-ML-method
 measuredSignal_ref (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 Measurement_package-class
 measuredSignal_ref, QuantitationType_assnref-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 mediaCol (zzz-ARRCompatibility-class),
 mediaCol
 measuredSignal_ref, QuantitationTypes_assnreflist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 mediaCol, PhysicalArray-method
 (zzz-ARRCompatibility-class),
 mediaFileGUID
 measuredSignal_ref, SourcesQuantitationType_assnreflist-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 (zzz-ARRCompatibility-class),
 mediaFileGUID
 measuredSignal_ref, TargetQuantitationType_assnref-method (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 (zzz-ARRCompatibility-class),
 mediaFileName
 MeasuredSignal_ref-class (zzz-DTTCompatibility-class),
 (zzz-DTTCompatibility-class), 13
 (zzz-ARRCompatibility-class),
 measurement 12

- mediaFileName,PhysicalArray-method
(zzz-ARRCompatibility-class),
[12](#)
- mediaRow (zzz-ARRCompatibility-class),
[12](#)
- mediaRow,PhysicalArray-method
(zzz-ARRCompatibility-class),
[12](#)
- mediaType (zzz-ARRCompatibility-class),
[12](#)
- mediaType,PhysicalArray-method
(zzz-ARRCompatibility-class),
[12](#)
- members_assnreflist
(zzz-DTTCompatibility-class),
[13](#)
- members_assnreflist,SecurityGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- Members_assnreflist-class
(zzz-DTTCompatibility-class),
[13](#)
- midInitials
(zzz-DTTCompatibility-class),
[13](#)
- midInitials,Person-method
(zzz-DTTCompatibility-class),
[13](#)
- mismatchInformation
(zzz-DTTCompatibility-class),
[13](#)
- mismatchInformation,MismatchInformation_assnlist-method
(zzz-DTTCompatibility-class),
[13](#)
- MismatchInformation-class
(zzz-DTTCompatibility-class),
[13](#)
- mismatchInformation_assnlist
(zzz-DTTCompatibility-class),
[13](#)
- mismatchInformation_assnlist,CompositePosition-method
(zzz-DTTCompatibility-class),
[13](#)
- mismatchInformation_assnlist,FeatureInformation-method
(zzz-DTTCompatibility-class),
[13](#)
- mismatchInformation_assnlist,ReporterPosition-method
(zzz-DTTCompatibility-class),
[13](#)
- MismatchInformation_assnlist-class
(zzz-DTTCompatibility-class),
[13](#)
- model (zzz-DTTCompatibility-class), [13](#)
- model,Hardware-method
(zzz-DTTCompatibility-class),
[13](#)
- name (zzz-DTTCompatibility-class), [13](#)
- name,Array-method
(zzz-DTTCompatibility-class),
[13](#)
- name,Array_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name,ArrayAttribute-method
(zzz-ARRCompatibility-class),
[12](#)
- name,ArrayDesign-method
(zzz-DTTCompatibility-class),
[13](#)
- name,ArrayDesign_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name,ArrayGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- name,ArrayGroup_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name,ArrayManufacture-method
(zzz-DTTCompatibility-class),
[13](#)
- name,ArrayManufacture_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name,BioAssay_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name,BioAssayCreation-method
(zzz-DTTCompatibility-class),
[13](#)
- name,BioAssayData_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name,BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
[13](#)
- name,BioAssayDataCluster_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name,BioAssayDimension-method
(zzz-DTTCompatibility-class),
[13](#)
- name,BioAssayDimension_ref-method
(zzz-DTTCompatibility-class),
[13](#)

- [13](#)
- name, BioAssayMap-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioAssayMap_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioAssayTreatment-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioMaterial_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioSample-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioSample_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioSequence-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioSequence_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioSource-method
(zzz-DTTCompatibility-class),
[13](#)
- name, BioSource_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Channel-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Channel_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeCompositeMap_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeGroup_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeSequence-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeSequence_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
[13](#)
- name, CompositeSequenceDimension_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Compound-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Compound_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, ConfidenceIndicator_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Contact_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Database-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Database_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DerivedBioAssay_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DerivedBioAssayData_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DerivedSignal-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DerivedSignal_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DesignElement_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DesignElementDimension_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, DesignElementMap_ref-method

- (zzz-DTTCompatibility-class),
13
- name,Error-method
(zzz-DTTCompatibility-class),
13
- name,Error_ref-method
(zzz-DTTCompatibility-class),
13
- name,ExpectedValue-method
(zzz-DTTCompatibility-class),
13
- name,ExpectedValue_ref-method
(zzz-DTTCompatibility-class),
13
- name,Experiment-method
(zzz-DTTCompatibility-class),
13
- name,ExperimentalFactor-method
(zzz-DTTCompatibility-class),
13
- name,ExperimentalFactor_ref-method
(zzz-DTTCompatibility-class),
13
- name,FactorValue-method
(zzz-DTTCompatibility-class),
13
- name,FactorValue_ref-method
(zzz-DTTCompatibility-class),
13
- name,Failed-method
(zzz-DTTCompatibility-class),
13
- name,Failed_ref-method
(zzz-DTTCompatibility-class),
13
- name,Feature-method
(zzz-DTTCompatibility-class),
13
- name,Feature_ref-method
(zzz-DTTCompatibility-class),
13
- name,FeatureDimension-method
(zzz-DTTCompatibility-class),
13
- name,FeatureDimension_ref-method
(zzz-DTTCompatibility-class),
13
- name,FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- name,FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- name,FeatureGroup_ref-method
(zzz-DTTCompatibility-class),
13
- name,FeatureReporterMap-method
(zzz-DTTCompatibility-class),
13
- name,FeatureReporterMap_ref-method
(zzz-DTTCompatibility-class),
13
- name,Hardware-method
(zzz-DTTCompatibility-class),
13
- name,Hardware_ref-method
(zzz-DTTCompatibility-class),
13
- name,Hybridization-method
(zzz-DTTCompatibility-class),
13
- name,Image-method
(zzz-DTTCompatibility-class),
13
- name,Image_ref-method
(zzz-DTTCompatibility-class),
13
- name,ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
- name,LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- name,LabeledExtract_ref-method
(zzz-DTTCompatibility-class),
13
- name,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- name,MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- name,MeasuredBioAssay_ref-method
(zzz-DTTCompatibility-class),
13
- name,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- name,MeasuredBioAssayData_ref-method
(zzz-DTTCompatibility-class),
13
- name,MeasuredSignal-method
(zzz-DTTCompatibility-class),
13

- name, MeasuredSignal_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, NameValueType-method
(zzz-DTTCompatibility-class),
[13](#)
- name, NodeValue-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Organization-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Organization_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Parameter-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Parameter_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Person-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Person_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PhysicalArrayDesign_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PhysicalBioAssay_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PresentAbsent-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PresentAbsent_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Protocol-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Protocol_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PValue-method
(zzz-DTTCompatibility-class),
[13](#)
- name, PValue_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, QuantitationType_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
[13](#)
- name, QuantitationTypeDimension_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
[13](#)
- name, QuantitationTypeMap_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Ratio-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Ratio_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Reporter-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Reporter_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
[13](#)
- name, ReporterCompositeMap_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, ReporterDimension-method
(zzz-DTTCompatibility-class),
[13](#)
- name, ReporterDimension_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, ReporterGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- name, ReporterGroup_ref-method
(zzz-DTTCompatibility-class),
[13](#)
- name, Security-method
(zzz-DTTCompatibility-class),

- 13
- name, Security_ref-method
(zzz-DTTCompatibility-class),
13
- name, SecurityGroup-method
(zzz-DTTCompatibility-class),
13
- name, SecurityGroup_ref-method
(zzz-DTTCompatibility-class),
13
- name, Software-method
(zzz-DTTCompatibility-class),
13
- name, Software_ref-method
(zzz-DTTCompatibility-class),
13
- name, SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- name, SpecializedQuantitationType_ref-method
(zzz-DTTCompatibility-class),
13
- name, StandardQuantitationType_ref-method
(zzz-DTTCompatibility-class),
13
- name, Transformation-method
(zzz-DTTCompatibility-class),
13
- name, Treatment-method
(zzz-DTTCompatibility-class),
13
- name, UserAttribute-method
(zzz-ARRCompatibility-class),
12
- name, Zone-method
(zzz-DTTCompatibility-class),
13
- name, Zone_ref-method
(zzz-DTTCompatibility-class),
13
- names, AffxArray-method
(zzz-NetAffxCompatibility-class),
14
- names, NetAffxResource-method
(NetAffxResource-class), 5
- namespace (zzz-ARRCompatibility-class),
12
- namespace, UserAttribute-method
(zzz-ARRCompatibility-class),
12
- nameValueType
(zzz-DTTCompatibility-class),
13
- nameValueType, PropertySets_assnlist-method
(zzz-DTTCompatibility-class),
13
- nameValueType, QualityControlStatistics_assnlist-method
(zzz-DTTCompatibility-class),
13
- nameValueType, SummaryStatistics_assnlist-method
(zzz-DTTCompatibility-class),
13
- NameValueType-class
(zzz-DTTCompatibility-class),
13
- NetAffxCompatibility, 6, 14
- NetAffxCompatibility-class
(AffyCompatible-class), 3
- NetAffxResource, 3, 4, 5, 9, 10, 15
- NetAffxResource-accessors (accessors), 2
- NetAffxResource-class, 5
- newSequence
(zzz-DTTCompatibility-class),
13
- newSequence, MismatchInformation-method
(zzz-DTTCompatibility-class),
13
- node (zzz-DTTCompatibility-class), 13
- node, Nodes_assnlist-method
(zzz-DTTCompatibility-class),
13
- Node-class
(zzz-DTTCompatibility-class),
13
- nodeContents
(zzz-DTTCompatibility-class),
13
- nodeContents, NodeContents_assnlist-method
(zzz-DTTCompatibility-class),
13
- NodeContents-class
(zzz-DTTCompatibility-class),
13
- nodeContents_assnlist
(zzz-DTTCompatibility-class),
13
- nodeContents_assnlist, Node-method
(zzz-DTTCompatibility-class),
13
- NodeContents_assnlist-class
(zzz-DTTCompatibility-class),
13
- nodes_assnlist
(zzz-DTTCompatibility-class),
13

- 13
- nodes_assnlist, BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
13
- nodes_assnlist, Node-method
(zzz-DTTCompatibility-class),
13
- Nodes_assnlist-class
(zzz-DTTCompatibility-class),
13
- nodeValue (zzz-DTTCompatibility-class),
13
- nodeValue, NodeValue_assnlist-method
(zzz-DTTCompatibility-class),
13
- NodeValue-class
(zzz-DTTCompatibility-class),
13
- nodeValue_assnlist
(zzz-DTTCompatibility-class),
13
- nodeValue_assnlist, Node-method
(zzz-DTTCompatibility-class),
13
- NodeValue_assnlist-class
(zzz-DTTCompatibility-class),
13
- normalizationDescription_assn
(zzz-DTTCompatibility-class),
13
- normalizationDescription_assn, ExperimentDesign-method
(zzz-DTTCompatibility-class),
13
- NormalizationDescription_assn-class
(zzz-DTTCompatibility-class),
13
- numArrays (zzz-DTTCompatibility-class),
13
- numArrays, ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- numberOfFeatures
(zzz-DTTCompatibility-class),
13
- numberOfFeatures, ArrayDesign-method
(zzz-DTTCompatibility-class),
13
- numberOfFeatures, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- numFeaturesPerCol
(zzz-DTTCompatibility-class),
13
- numFeaturesPerCol, ZoneLayout-method
(zzz-DTTCompatibility-class),
13
- numFeaturesPerRow
(zzz-DTTCompatibility-class),
13
- numFeaturesPerRow, ZoneLayout-method
(zzz-DTTCompatibility-class),
13
- ontologyEntries_assnlist
(zzz-DTTCompatibility-class),
13
- ontologyEntries_assnlist, BioSequence-method
(zzz-DTTCompatibility-class),
13
- OntologyEntries_assnlist-class
(zzz-DTTCompatibility-class),
13
- ontologyEntry
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Action_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Annotations_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Associations_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Category_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Characteristics_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, CompoundIndices_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, ControlType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, DataType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, DefectType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, FailTypes_assnlist-method
(zzz-DTTCompatibility-class),
13

- ontologyEntry, FeatureShape_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, FiducialType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Format_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, MaterialType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, OntologyEntries_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Parameters_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, PolymerType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Roles_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Scale_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Species_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, SubstrateType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, SurfaceType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, TechnologyType_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Type_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Types_assnlist-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, Value_assn-method
(zzz-DTTCompatibility-class),
13
- ontologyEntry, WarningType_assn-method
(zzz-DTTCompatibility-class),
13
- OntologyEntry-class
(zzz-DTTCompatibility-class),
13
- ontologyReference_assn
(zzz-DTTCompatibility-class),
13
- ontologyReference_assn, OntologyEntry-method
(zzz-DTTCompatibility-class),
13
- OntologyReference_assn-class
(zzz-DTTCompatibility-class),
13
- organization
(zzz-DTTCompatibility-class),
13
- organization, Contact_assnlist-method
(zzz-DTTCompatibility-class),
13
- organization, Organization_assnlist-method
(zzz-DTTCompatibility-class),
13
- Organization-class
(zzz-DTTCompatibility-class),
13
- Organization_assnlist-class
(zzz-DTTCompatibility-class),
13
- organization_ref
(zzz-DTTCompatibility-class),
13
- organization_ref, Affiliation_assnref-method
(zzz-DTTCompatibility-class),
13
- organization_ref, ArrayManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, Contacts_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, DesignProviders_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, HardwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, Members_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, Owner_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, Parent_assnref-method
(zzz-DTTCompatibility-class),

- 13
- organization_ref, Performer_assnref-method
(zzz-DTTCompatibility-class),
13
- organization_ref, Providers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, SoftwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- organization_ref, SourceContact_assnreflist-method
(zzz-DTTCompatibility-class),
13
- Organization_ref-class
(zzz-DTTCompatibility-class),
13
- orientationMark
(zzz-DTTCompatibility-class),
13
- orientationMark, ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- orientationMarkPosition
(zzz-DTTCompatibility-class),
13
- orientationMarkPosition, ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- originalProjectName
(zzz-ARRCompatibility-class),
12
- originalProjectName, ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- originRelativeTo
(zzz-DTTCompatibility-class),
13
- originRelativeTo, Array-method
(zzz-DTTCompatibility-class),
13
- otherKind (zzz-DTTCompatibility-class),
13
- otherKind, Measurement-method
(zzz-DTTCompatibility-class),
13
- owner_assnreflist
(zzz-DTTCompatibility-class),
13
- owner_assnreflist, Security-method
(zzz-DTTCompatibility-class),
13
- Owner_assnreflist-class
(zzz-DTTCompatibility-class),
13
- pages (zzz-DTTCompatibility-class), 13
- pages, BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- parameter (zzz-DTTCompatibility-class),
13
- parameter, ParameterTypes_assnlist-method
(zzz-DTTCompatibility-class),
13
- Parameter-class
(zzz-DTTCompatibility-class),
13
- parameter_ref
(zzz-DTTCompatibility-class),
13
- parameter_ref, ParameterType_assnref-method
(zzz-DTTCompatibility-class),
13
- Parameter_ref-class
(zzz-DTTCompatibility-class),
13
- parameters_assnlist
(zzz-DTTCompatibility-class),
13
- parameters_assnlist, BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- Parameters_assnlist-class
(zzz-DTTCompatibility-class),
13
- parameterType_assnref
(zzz-DTTCompatibility-class),
13
- parameterType_assnref, ParameterValue-method
(zzz-DTTCompatibility-class),
13
- ParameterType_assnref-class
(zzz-DTTCompatibility-class),
13
- parameterTypes_assnlist
(zzz-DTTCompatibility-class),
13
- parameterTypes_assnlist, Hardware-method
(zzz-DTTCompatibility-class),
13
- parameterTypes_assnlist, Protocol-method
(zzz-DTTCompatibility-class),
13
- parameterTypes_assnlist, Software-method
(zzz-DTTCompatibility-class),

- 13
- ParameterTypes_assnlist-class
(zzz-DTTCompatibility-class),
13
- parameterValue
(zzz-DTTCompatibility-class),
13
- parameterValue,ParameterValues_assnlist-method
(zzz-DTTCompatibility-class),
13
- ParameterValue-class
(zzz-DTTCompatibility-class),
13
- parameterValues_assnlist
(zzz-DTTCompatibility-class),
13
- parameterValues_assnlist,HardwareApplication-method
(zzz-DTTCompatibility-class),
13
- parameterValues_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- parameterValues_assnlist,SoftwareApplication-method
(zzz-DTTCompatibility-class),
13
- ParameterValues_assnlist-class
(zzz-DTTCompatibility-class),
13
- parent_assnref
(zzz-DTTCompatibility-class),
13
- parent_assnref,Organization-method
(zzz-DTTCompatibility-class),
13
- Parent_assnref-class
(zzz-DTTCompatibility-class),
13
- password (accessors), 2
- password,NetAffxResource-method
(NetAffxResource-class), 5
- patassignmentMethod
(zzz-ARRCompatibility-class),
12
- patassignmentMethod,PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- performer_assnref
(zzz-DTTCompatibility-class),
13
- performer_assnref,Audit-method
(zzz-DTTCompatibility-class),
13
- Performer_assnref-class
(zzz-DTTCompatibility-class),
13
- performers_assnreflist
(zzz-DTTCompatibility-class),
13
- performers_assnreflist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- Performers_assnreflist-class
(zzz-DTTCompatibility-class),
13
- Person-class
(zzz-DTTCompatibility-class),
13
- Person_assnlist-class
(zzz-DTTCompatibility-class),
13
- person_ref
(zzz-DTTCompatibility-class),
13
- person_ref,ArrayManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,Contacts_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,DesignProviders_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,HardwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,Members_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,Owner_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,Performer_assnref-method
(zzz-DTTCompatibility-class),
13
- person_ref,Performers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,Providers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,SoftwareManufacturers_assnreflist-method
(zzz-DTTCompatibility-class),
13
- person_ref,SourceContact_assnreflist-method

- (zzz-DTTCompatibility-class),
13
- Person_ref-class
(zzz-DTTCompatibility-class),
13
- phone (zzz-DTTCompatibility-class), 13
- phone, Organization-method
(zzz-DTTCompatibility-class),
13
- phone, Person-method
(zzz-DTTCompatibility-class),
13
- physicalArray
(zzz-ARRCompatibility-class),
12
- physicalArray, PhysicalArrays-method
(zzz-ARRCompatibility-class),
12
- PhysicalArray-class
(zzz-ARRCompatibility-class),
12
- physicalArrayDesign
(zzz-DTTCompatibility-class),
13
- physicalArrayDesign, ArrayDesign_assnlist-method
(zzz-DTTCompatibility-class),
13
- PhysicalArrayDesign-class
(zzz-DTTCompatibility-class),
13
- physicalArrayDesign_ref
(zzz-DTTCompatibility-class),
13
- physicalArrayDesign_ref, ArrayDesign_assnref-method
(zzz-DTTCompatibility-class),
13
- PhysicalArrayDesign_ref-class
(zzz-DTTCompatibility-class),
13
- physicalArrays
(zzz-ARRCompatibility-class),
12
- physicalArrays, ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- PhysicalArrays-class
(zzz-ARRCompatibility-class),
12
- physicalBioAssay
(zzz-DTTCompatibility-class),
13
- physicalBioAssay, BioAssay_assnlist-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay, PhysicalBioAssay_assnlist-method
(zzz-DTTCompatibility-class),
13
- PhysicalBioAssay-class
(zzz-DTTCompatibility-class),
13
- PhysicalBioAssay_assnlist-class
(zzz-DTTCompatibility-class),
13
- PhysicalBioAssay_assnref-class
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, BioAssay_assnref-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, BioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, PhysicalBioAssay_assnref-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, PhysicalBioAssaySource_assnref-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, PhysicalBioAssayTarget_assnref-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, SourceBioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, Target_assnref-method
(zzz-DTTCompatibility-class),
13
- physicalBioAssay_ref, TopLevelBioAssays_assnreflist-method
(zzz-DTTCompatibility-class),
13
- PhysicalBioAssay_ref-class
(zzz-DTTCompatibility-class),
13
- physicalBioAssayData_assnlist
(zzz-DTTCompatibility-class),
13
- physicalBioAssayData_assnlist, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- PhysicalBioAssayData_assnlist-class
(zzz-DTTCompatibility-class),

- 13
- physicalBioAssaySource_assnref
(zzz-DTTCompatibility-class),
13
- physicalBioAssaySource_assnref,FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- PhysicalBioAssaySource_assnref-class
(zzz-DTTCompatibility-class),
13
- PhysicalBioAssayTarget_assnref-class
(zzz-DTTCompatibility-class),
13
- polymerType_assn
(zzz-DTTCompatibility-class),
13
- polymerType_assn,BioSequence-method
(zzz-DTTCompatibility-class),
13
- PolymerType_assn-class
(zzz-DTTCompatibility-class),
13
- position (zzz-DTTCompatibility-class),
13
- position,Position_assn-method
(zzz-DTTCompatibility-class),
13
- Position-class
(zzz-DTTCompatibility-class),
13
- position_assn
(zzz-DTTCompatibility-class),
13
- position_assn,Feature-method
(zzz-DTTCompatibility-class),
13
- position_assn,Fiducial-method
(zzz-DTTCompatibility-class),
13
- Position_assn-class
(zzz-DTTCompatibility-class),
13
- positionDelta
(zzz-DTTCompatibility-class),
13
- positionDelta,PositionDelta_assn-method
(zzz-DTTCompatibility-class),
13
- PositionDelta-class
(zzz-DTTCompatibility-class),
13
- positionDelta_assn
(zzz-DTTCompatibility-class),
13
- positionDelta_assn,FeatureDefect-method
(zzz-DTTCompatibility-class),
13
- positionDelta_assn,ZoneDefect-method
(zzz-DTTCompatibility-class),
13
- PositionDelta_assn-class
(zzz-DTTCompatibility-class),
13
- presentAbsent
(zzz-DTTCompatibility-class),
13
- presentAbsent,QuantitationType_assnlist-method
(zzz-DTTCompatibility-class),
13
- PresentAbsent-class
(zzz-DTTCompatibility-class),
13
- presentAbsent_ref
(zzz-DTTCompatibility-class),
13
- presentAbsent_ref,QuantitationType_assnref-method
(zzz-DTTCompatibility-class),
13
- presentAbsent_ref,QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class),
13
- presentAbsent_ref,SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class),
13
- presentAbsent_ref,TargetQuantitationType_assnref-method
(zzz-DTTCompatibility-class),
13
- PresentAbsent_ref-class
(zzz-DTTCompatibility-class),
13
- producerTransformation_assn
(zzz-DTTCompatibility-class),
13
- producerTransformation_assn,DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
13
- ProducerTransformation_assn-class
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Array-method
(zzz-DTTCompatibility-class),

- 13
propertySets_assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,ArrayGroup-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,ArrayManufacture-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,ArrayManufactureDeviation-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,Audit-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BibliographicReference-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioAssayCreation-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioAssayDimension-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioAssayMap-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioAssayMapping-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioAssayTreatment-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioDataCube-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioDataTuples-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioMaterialMeasurement-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioSample-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioSequence-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,BioSource-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,Channel-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,CompositeGroup-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,CompositePosition-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,CompositeSequence-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,Compound-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,CompoundMeasurement-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,ConcentrationUnit-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,Database-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,DatabaseEntry-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,DerivedBioAssay-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,DerivedSignal-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,Description-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,DesignElementMapping-method
(zzz-DTTCompatibility-class),
13
propertySets_assnlist,DistanceUnit-method

- (zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Error-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ExpectedValue-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Experiment-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ExperimentalFactor-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ExperimentDesign-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ExternalReference-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FactorValue-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Failed-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Feature-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FeatureDefect-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FeatureDimension-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FeatureInformation-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FeatureLocation-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Fiducial-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Hardware-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,HardwareApplication-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Hybridization-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Image-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ManufactureLIMS-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,MassUnit-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Measurement-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,MismatchInformation-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,NameValueType-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Node-method
(zzz-DTTCompatibility-class),
13

- propertySets_assnlist,NodeContents-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,NodeValue-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,OntologyEntry-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Organization-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Parameter-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ParameterValue-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Person-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Position-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,PositionDelta-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Protocol-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,PValue-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,QuantitationTypeMapping-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,QuantityUnit-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Ratio-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Reporter-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ReporterDimension-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ReporterPosition-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Security-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,SecurityGroup-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,SeqFeature-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,SeqFeatureLocation-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,SequencePosition-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Software-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,SoftwareApplication-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,TemperatureUnit-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,TimeUnit-method
(zzz-DTTCompatibility-class),
13

- 13
- propertySets_assnlist,Transformation-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Treatment-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,VolumeUnit-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,Zone-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ZoneDefect-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ZoneGroup-method
(zzz-DTTCompatibility-class),
13
- propertySets_assnlist,ZoneLayout-method
(zzz-DTTCompatibility-class),
13
- PropertySets_assnlist-class
(zzz-DTTCompatibility-class),
13
- protocol (zzz-DTTCompatibility-class),
13
- protocol,Protocol_assnlist-method
(zzz-DTTCompatibility-class),
13
- Protocol-class
(zzz-DTTCompatibility-class),
13
- protocol_assnlist
(zzz-DTTCompatibility-class),
13
- protocol_assnlist,Protocol_package-method
(zzz-DTTCompatibility-class),
13
- Protocol_assnlist-class
(zzz-DTTCompatibility-class),
13
- protocol_assnref
(zzz-DTTCompatibility-class),
13
- protocol_assnref,ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- Protocol_assnref-class
(zzz-DTTCompatibility-class),
13
- protocol_package
(zzz-DTTCompatibility-class),
13
- protocol_package,MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- Protocol_package-class
(zzz-DTTCompatibility-class),
13
- protocol_ref
(zzz-DTTCompatibility-class),
13
- protocol_ref,Protocol_assnref-method
(zzz-DTTCompatibility-class),
13
- Protocol_ref-class
(zzz-DTTCompatibility-class),
13
- protocolApplication
(zzz-DTTCompatibility-class),
13
- protocolApplication,ProtocolApplications_assnlist-method
(zzz-DTTCompatibility-class),
13
- ProtocolApplication-class
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,ArrayManufacture-method
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,BioAssayCreation-method
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,BioAssayMap-method
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,BioAssayTreatment-method
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- protocolApplications_assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class),

- 13 PValue-class
- protocolApplications_assnlist,Hybridization-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13 pvalue_ref
- protocolApplications_assnlist,ImageAcquisition-method(zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13 pvalue_ref,ConfidenceIndicators_assnreflist-method
- protocolApplications_assnlist,PhysicalArrayDesign-method(zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13 pvalue_ref,QuantitationType_assnref-method
- protocolApplications_assnlist,QuantitationTypeMap-method(zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13 pvalue_ref,QuantitationTypes_assnreflist-method
- protocolApplications_assnlist,ReporterCompositeMap-method(zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13 pvalue_ref,SourcesQuantitationType_assnreflist-method
- protocolApplications_assnlist,Transformation-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13 pvalue_ref,TargetQuantitationType_assnref-method
- protocolApplications_assnlist,Treatment-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13 PValue_ref-class
- ProtocolApplications_assnlist-class (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- 13
- providers_assnreflist quality (zzz-DTTCompatibility-class), 13
- (zzz-DTTCompatibility-class), quality,ManufactureLIMS-method
13 (zzz-DTTCompatibility-class),
13
- providers_assnreflist,Experiment-method quality,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), (zzz-DTTCompatibility-class),
13 13
- Providers_assnreflist-class qualityControlDescription_assn
(zzz-DTTCompatibility-class), (zzz-DTTCompatibility-class),
13 13
- publication qualityControlDescription_assn,ExperimentDesign-method
(zzz-DTTCompatibility-class), (zzz-DTTCompatibility-class),
13 13
- publication,BibliographicReference-method QualityControlDescription_assn-class
(zzz-DTTCompatibility-class), (zzz-DTTCompatibility-class),
13 13
- publisher (zzz-DTTCompatibility-class), qualityControlStatistics_assnlist
13 (zzz-DTTCompatibility-class),
13
- publisher,BibliographicReference-method qualityControlStatistics_assnlist,ArrayManufacture-meth
(zzz-DTTCompatibility-class), (zzz-DTTCompatibility-class),
13 13
- pvalue (zzz-DTTCompatibility-class), 13
- pvalue,ConfidenceIndicator_assnlist-method qualityControlStatistics_assnlist,BioSample-method
(zzz-DTTCompatibility-class), (zzz-DTTCompatibility-class),
13 13
- pvalue,QuantitationType_assnlist-method qualityControlStatistics_assnlist,BioSource-method
(zzz-DTTCompatibility-class), (zzz-DTTCompatibility-class),
13 13

- qualityControlStatistics_assnlist, LabeledExtract-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QualityControlStatistics_assnlist-class
(zzz-DTTCompatibility-class), 13
- quantitationDimension_assnref
(zzz-DTTCompatibility-class), 13
- quantitationDimension_assnref, NodeContents-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationDimension_assnref-class
(zzz-DTTCompatibility-class), 13
- quantitationType_assnlist
(zzz-DTTCompatibility-class), 13
- quantitationType_assnlist, QuantitationType_package-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationType_assnlist-class
(zzz-DTTCompatibility-class), 13
- quantitationType_assnref
(zzz-DTTCompatibility-class), 13
- quantitationType_assnref, QuantitationTypeTuple-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationType_assnref-class
(zzz-DTTCompatibility-class), 13
- quantitationType_package
(zzz-DTTCompatibility-class), 13
- quantitationType_package, MAGe-ML-method
(zzz-DTTCompatibility-class), 13
- QuantitationType_package-class
(zzz-DTTCompatibility-class), 13
- QuantitationType_ref-class
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension, QuantitationTypeDimension_assnlist-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationTypeDimension-class
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_assnlist
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_assnlist, BioAssayData_package-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationTypeDimension_assnlist-class
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_assnref
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_assnref, DerivedBioAssayData-package-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_assnref, MeasuredBioAssayData-package-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationTypeDimension_assnref-class
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_ref
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_ref, QuantitationDimension_assnlist-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_ref, QuantitationTypeDimension_assnlist-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationTypeDimension_ref-class
(zzz-DTTCompatibility-class), 13
- quantitationTypeMap
(zzz-DTTCompatibility-class), 13
- quantitationTypeMap, QuantitationTypeMap_assnlist-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationTypeMap-class
(zzz-DTTCompatibility-class), 13
- quantitationTypeMap_assnlist
(zzz-DTTCompatibility-class), 13
- quantitationTypeMap_assnlist, BioAssayData_package-method (zzz-DTTCompatibility-class),
(zzz-DTTCompatibility-class), 13
- QuantitationTypeMap_assnlist-class
(zzz-DTTCompatibility-class), 13

[13](#)
 quantitationTypeMap_ref
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMap_ref, QuantitationTypeMaps_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantitationTypeMap_ref-class
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMapping
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMapping, QuantitationTypeMapping_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantitationTypeMapping-class
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMapping_assn
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMapping_assn, Transformation-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantitationTypeMapping_assn-class
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, DerivedSignal-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, Error-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, ExpectedValue-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, Failed-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, MeasuredSignal-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, PresentAbsent-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, PValue-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, QuantitationTypeMapping-assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, Ratio-method
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeMaps_assnreflist, SpecializedQuantitationTypeMapping-assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantitationTypeMaps_assnreflist-class
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypes_assnreflist
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypes_assnreflist, QuantitationTypeDimensionality-assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantitationTypes_assnreflist-class
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeTuple
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeTuple, QuantitationTypeTuples_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantitationTypeTuple-class
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeTuples_assnreflist
 (zzz-DTTCompatibility-class),
[13](#)
 quantitationTypeTuples_assnreflist, DesignElementTuple-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantitationTypeTuples_assnreflist-class
 (zzz-DTTCompatibility-class),
[13](#)
 quantityUnit
 (zzz-DTTCompatibility-class),
[13](#)
 quantityUnit, Unit_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)
 QuantityUnit-class
 (zzz-DTTCompatibility-class),
[13](#)
 ratio (zzz-DTTCompatibility-class), [13](#)
 ratio, QuantitationType_assnreflist-method
 (zzz-DTTCompatibility-class),
[13](#)

- Ratio-class (zzz-DTTCompatibility-class), 13
- ratio_ref (zzz-DTTCompatibility-class), 13
- ratio_ref, QuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- ratio_ref, QuantitationTypes_assnreflist-method (zzz-DTTCompatibility-class), 13
- ratio_ref, SourcesQuantitationType_assnreflist-method (zzz-DTTCompatibility-class), 13
- ratio_ref, TargetQuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- Ratio_ref-class (zzz-DTTCompatibility-class), 13
- readAnnotation (NetAffxResource-class), 5
- readAnnotation, NetAffxResource, character, character-methods (NetAffxResource-class), 5
- readAnnotation, NetAffxResource, character, numeric-methods (NetAffxResource-class), 5
- readAnnotation, NetAffxResource, missing, AffxAnnotation-method (NetAffxResource-class), 5
- readAnnotation, NetAffxResource, numeric, character-methods (NetAffxResource-class), 5
- readAnnotation, NetAffxResource, numeric, numeric-methods (NetAffxResource-class), 5
- readArr, 8
- readMage, 8
- readNetAffx, 9
- readXml, 8, 9, 10, 11
- regions_assnlist (zzz-DTTCompatibility-class), 13
- regions_assnlist, SeqFeature-method (zzz-DTTCompatibility-class), 13
- Regions_assnlist-class (zzz-DTTCompatibility-class), 13
- releaseDate (zzz-DTTCompatibility-class), 13
- releaseDate, SoftwareApplication-method (zzz-DTTCompatibility-class), 13
- replacedLength (zzz-DTTCompatibility-class), 13
- replacedLength, MismatchInformation-method (zzz-DTTCompatibility-class), 13
- replicateDescription_assn (zzz-DTTCompatibility-class), 13
- replicateDescription_assn, ExperimentDesign-method (zzz-DTTCompatibility-class), 13
- Reporter-class (zzz-DTTCompatibility-class), 13
- Reporter-methods (zzz-DTTCompatibility-class), 13
- Reporter_assnlist (zzz-DTTCompatibility-class), 13
- Reporter_assnlist, DesignElement_assnlist-method (zzz-DTTCompatibility-class), 13
- Reporter_assnlist, Reporter_assnlist-method (zzz-DTTCompatibility-class), 13
- Reporter_assnlist-class (zzz-DTTCompatibility-class), 13
- Reporter_assnlist, DesignElement_package-method (zzz-DTTCompatibility-class), 13
- Reporter_assnlist-method (zzz-DTTCompatibility-class), 13
- Reporter_assnlist, ReporterPosition-method (zzz-DTTCompatibility-class), 13
- Reporter_assnlist-class (zzz-DTTCompatibility-class), 13
- Reporter_assnref (zzz-DTTCompatibility-class), 13
- Reporter_assnref, FeatureReporterMap-method (zzz-DTTCompatibility-class), 13
- Reporter_assnref, ReporterPosition-method (zzz-DTTCompatibility-class), 13
- Reporter_assnref-class (zzz-DTTCompatibility-class), 13
- Reporter_ref (zzz-DTTCompatibility-class), 13
- Reporter_ref, DesignElement_assnref-method (zzz-DTTCompatibility-class), 13

- reporter_ref,Reporter_assnref-method (zzz-DTTCompatibility-class), 13
- reporter_ref,Reporters_assnreflist-method (zzz-DTTCompatibility-class), 13
- Reporter_ref-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMap (zzz-DTTCompatibility-class), 13
- reporterCompositeMap,DesignElementMap_assnlist-method(zzz-DTTCompatibility-class), 13
- reporterCompositeMap,ReporterCompositeMap_assnlist-method(zzz-DTTCompatibility-class), 13
- ReporterCompositeMap-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_assnlist (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_assnlist,DesignElement_package-method(zzz-DTTCompatibility-class), 13
- ReporterCompositeMap_assnlist-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_ref (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_ref,DesignElementMaps_assnreflist(zzz-DTTCompatibility-class), 13
- reporterCompositeMap_ref,ReporterCompositeMaps_assnreflist(zzz-DTTCompatibility-class), 13
- ReporterCompositeMap_ref-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMaps_assnreflist (zzz-DTTCompatibility-class), 13
- reporterCompositeMaps_assnreflist,CompositeSequence-method(zzz-DTTCompatibility-class), 13
- ReporterCompositeMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- reporterDimension (zzz-DTTCompatibility-class), 13
- reporterDimension,DesignElementDimension_assnlist-method(zzz-DTTCompatibility-class), 13
- ReporterDimension-class (zzz-DTTCompatibility-class), 13
- reporterDimension_ref (zzz-DTTCompatibility-class), 13
- reporterDimension_ref,DesignElementDimension_assnref-method(zzz-DTTCompatibility-class), 13
- ReporterDimension_ref-class (zzz-DTTCompatibility-class), 13
- reporterGroup (zzz-DTTCompatibility-class), 13
- reporterGroup,ReporterGroup_assnlist-method (zzz-DTTCompatibility-class), 13
- ReporterGroup-class (zzz-DTTCompatibility-class), 13
- reporterGroup_assnlist (zzz-DTTCompatibility-class), 13
- reporterGroup_assnlist,ArrayDesign_package-method (zzz-DTTCompatibility-class), 13
- ReporterGroup_assnlist-class (zzz-DTTCompatibility-class), 13
- reporterGroup_ref (zzz-DTTCompatibility-class), 13
- reporterGroup_ref,ReporterGroups_assnreflist-method (zzz-DTTCompatibility-class), 13
- ReporterGroup_ref-class (zzz-DTTCompatibility-class), 13
- reporterGroups_assnreflist (zzz-DTTCompatibility-class), 13
- reporterGroups_assnreflist,ArrayDesign-method (zzz-DTTCompatibility-class), 13
- reporterGroups_assnreflist,PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 13

- 13
- ReporterGroups_assnreflist-class
(zzz-DTTCompatibility-class),
13
- reporterPosition
(zzz-DTTCompatibility-class),
13
- reporterPosition,Coordinate_assn-method
(zzz-DTTCompatibility-class),
13
- reporterPosition,ReporterPositionSources_assnlist-method
(zzz-DTTCompatibility-class),
13
- ReporterPosition-class
(zzz-DTTCompatibility-class),
13
- reporterPositionSources_assnlist
(zzz-DTTCompatibility-class),
13
- reporterPositionSources_assnlist,ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
13
- ReporterPositionSources_assnlist-class
(zzz-DTTCompatibility-class),
13
- reporters_assnreflist
(zzz-DTTCompatibility-class),
13
- reporters_assnreflist,ReporterDimension-method
(zzz-DTTCompatibility-class),
13
- reporters_assnreflist,ReporterGroup-method
(zzz-DTTCompatibility-class),
13
- Reporters_assnreflist-class
(zzz-DTTCompatibility-class),
13
- required(zzz-ARRCompatibility-class),
12
- required,UserAttribute-method
(zzz-ARRCompatibility-class),
12
- roles_assnlist
(zzz-DTTCompatibility-class),
13
- roles_assnlist,Organization-method
(zzz-DTTCompatibility-class),
13
- roles_assnlist,Person-method
(zzz-DTTCompatibility-class),
13
- Roles_assnlist-class
(zzz-DTTCompatibility-class),
13
- scale_assn
(zzz-DTTCompatibility-class),
13
- scale_assn,DerivedSignal-method
(zzz-DTTCompatibility-class),
13
- scale_assn,Error-method
(zzz-DTTCompatibility-class),
13
- scale_assn,ExpectedValue-method
(zzz-DTTCompatibility-class),
13
- scale_assn,Failed-method
(zzz-DTTCompatibility-class),
13
- scale_assn,MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- scale_assn,NodeValue-method
(zzz-DTTCompatibility-class),
13
- scale_assn,PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- scale_assn,PValue-method
(zzz-DTTCompatibility-class),
13
- scale_assn,Ratio-method
(zzz-DTTCompatibility-class),
13
- scale_assn,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- Scale_assn-class
(zzz-DTTCompatibility-class),
13
- security(zzz-DTTCompatibility-class),
13
- security,Security_assnlist-method
(zzz-DTTCompatibility-class),
13
- Security-class
(zzz-DTTCompatibility-class),
13
- security_assnlist
(zzz-DTTCompatibility-class),
13
- security_assnlist,AuditAndSecurity_package-method
(zzz-DTTCompatibility-class),
13

- Security_assnlist-class
(zzz-DTTCompatibility-class),
13
- security_assnref
(zzz-DTTCompatibility-class),
13
- security_assnref,Array-method
(zzz-DTTCompatibility-class),
13
- security_assnref,ArrayDesign-method
(zzz-DTTCompatibility-class),
13
- security_assnref,ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- security_assnref,ArrayManufacture-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Audit-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioAssayCreation-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioAssayDataCluster-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioAssayDimension-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioAssayMap-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioAssayTreatment-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioSample-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioSequence-method
(zzz-DTTCompatibility-class),
13
- security_assnref,BioSource-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Channel-method
(zzz-DTTCompatibility-class),
13
- security_assnref,CompositeCompositeMap-method
(zzz-DTTCompatibility-class),
13
- security_assnref,CompositeGroup-method
(zzz-DTTCompatibility-class),
13
- security_assnref,CompositeSequence-method
(zzz-DTTCompatibility-class),
13
- security_assnref,CompositeSequenceDimension-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Compound-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Database-method
(zzz-DTTCompatibility-class),
13
- security_assnref,DerivedBioAssay-method
(zzz-DTTCompatibility-class),
13
- security_assnref,DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
13
- security_assnref,DerivedSignal-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Description-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Error-method
(zzz-DTTCompatibility-class),
13
- security_assnref,ExpectedValue-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Experiment-method
(zzz-DTTCompatibility-class),
13
- security_assnref,ExperimentalFactor-method
(zzz-DTTCompatibility-class),
13
- security_assnref,ExperimentDesign-method
(zzz-DTTCompatibility-class),
13
- security_assnref,FactorValue-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Failed-method
(zzz-DTTCompatibility-class),
13
- security_assnref,Feature-method
(zzz-DTTCompatibility-class),

- 13
- security_assnref, FeatureDimension-method
(zzz-DTTCompatibility-class),
13
- security_assnref, FeatureExtraction-method
(zzz-DTTCompatibility-class),
13
- security_assnref, FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- security_assnref, FeatureReporterMap-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Fiducial-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Hardware-method
(zzz-DTTCompatibility-class),
13
- security_assnref, HardwareApplication-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Hybridization-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Image-method
(zzz-DTTCompatibility-class),
13
- security_assnref, ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
- security_assnref, LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- security_assnref, MAGE-ML-method
(zzz-DTTCompatibility-class),
13
- security_assnref, ManufactureLIMS-method
(zzz-DTTCompatibility-class),
13
- security_assnref, ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class),
13
- security_assnref, MeasuredBioAssay-method
(zzz-DTTCompatibility-class),
13
- security_assnref, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- security_assnref, MeasuredSignal-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Node-method
(zzz-DTTCompatibility-class),
13
- security_assnref, NodeContents-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Organization-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Parameter-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Person-method
(zzz-DTTCompatibility-class),
13
- security_assnref, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- security_assnref, PhysicalBioAssay-method
(zzz-DTTCompatibility-class),
13
- security_assnref, PresentAbsent-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Protocol-method
(zzz-DTTCompatibility-class),
13
- security_assnref, ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- security_assnref, PValue-method
(zzz-DTTCompatibility-class),
13
- security_assnref, QuantitationTypeDimension-method
(zzz-DTTCompatibility-class),
13
- security_assnref, QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Ratio-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Reporter-method
(zzz-DTTCompatibility-class),
13
- security_assnref, ReporterCompositeMap-method
(zzz-DTTCompatibility-class),
13
- security_assnref, ReporterDimension-method
(zzz-DTTCompatibility-class),
13
- security_assnref, ReporterGroup-method

- (zzz-DTTCompatibility-class),
13
- security_assnref, Security-method
(zzz-DTTCompatibility-class),
13
- security_assnref, SecurityGroup-method
(zzz-DTTCompatibility-class),
13
- security_assnref, SeqFeature-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Software-method
(zzz-DTTCompatibility-class),
13
- security_assnref, SoftwareApplication-method
(zzz-DTTCompatibility-class),
13
- security_assnref, SpecializedQuantitationType-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Transformation-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Treatment-method
(zzz-DTTCompatibility-class),
13
- security_assnref, Zone-method
(zzz-DTTCompatibility-class),
13
- Security_assnref-class
(zzz-DTTCompatibility-class),
13
- security_ref
(zzz-DTTCompatibility-class),
13
- security_ref, Security_assnref-method
(zzz-DTTCompatibility-class),
13
- Security_ref-class
(zzz-DTTCompatibility-class),
13
- securityGroup
(zzz-DTTCompatibility-class),
13
- securityGroup, SecurityGroup_assnlist-method
(zzz-DTTCompatibility-class),
13
- SecurityGroup-class
(zzz-DTTCompatibility-class),
13
- securityGroup_assnlist
(zzz-DTTCompatibility-class),
13
- securityGroup_assnlist, AuditAndSecurity_package-method
(zzz-DTTCompatibility-class),
13
- SecurityGroup_assnlist-class
(zzz-DTTCompatibility-class),
13
- securityGroup_ref
(zzz-DTTCompatibility-class),
13
- securityGroup_ref, SecurityGroups_assnreflist-method
(zzz-DTTCompatibility-class),
13
- SecurityGroup_ref-class
(zzz-DTTCompatibility-class),
13
- securityGroups_assnreflist
(zzz-DTTCompatibility-class),
13
- securityGroups_assnreflist, Security-method
(zzz-DTTCompatibility-class),
13
- SecurityGroups_assnreflist-class
(zzz-DTTCompatibility-class),
13
- seqFeature
(zzz-DTTCompatibility-class),
13
- seqFeature, SeqFeatures_assnlist-method
(zzz-DTTCompatibility-class),
13
- SeqFeature-class
(zzz-DTTCompatibility-class),
13
- seqFeatureLocation
(zzz-DTTCompatibility-class),
13
- seqFeatureLocation, Regions_assnlist-method
(zzz-DTTCompatibility-class),
13
- seqFeatureLocation, Subregions_assnlist-method
(zzz-DTTCompatibility-class),
13
- SeqFeatureLocation-class
(zzz-DTTCompatibility-class),
13
- seqFeatures_assnlist
(zzz-DTTCompatibility-class),
13
- seqFeatures_assnlist, BioSequence-method
(zzz-DTTCompatibility-class),
13

- SeqFeatures_assnlist-class
(zzz-DTTCompatibility-class),
13
- sequenceDatabases_assnlist
(zzz-DTTCompatibility-class),
13
- sequenceDatabases_assnlist,BioSequence-method
(zzz-DTTCompatibility-class),
13
- SequenceDatabases_assnlist-class
(zzz-DTTCompatibility-class),
13
- sequencePosition
(zzz-DTTCompatibility-class),
13
- sequencePosition,Coordinate_assn-method
(zzz-DTTCompatibility-class),
13
- SequencePosition-class
(zzz-DTTCompatibility-class),
13
- serialNumber
(zzz-DTTCompatibility-class),
13
- serialNumber,HardwareApplication-method
(zzz-DTTCompatibility-class),
13
- show,.TypedAssociation-method
(Internal-class), 4
- show,AffyCompatible-method
(AffyCompatible-class), 3
- show,NetAffxResource-method
(NetAffxResource-class), 5
- software (zzz-DTTCompatibility-class),
13
- software,Software_assnlist-method
(zzz-DTTCompatibility-class),
13
- Software-class
(zzz-DTTCompatibility-class),
13
- software_assnlist
(zzz-DTTCompatibility-class),
13
- software_assnlist,Protocol_package-method
(zzz-DTTCompatibility-class),
13
- Software_assnlist-class
(zzz-DTTCompatibility-class),
13
- software_assnref
(zzz-DTTCompatibility-class),
13
- software_assnref,SoftwareApplication-method
(zzz-DTTCompatibility-class),
13
- Software_assnref-class
(zzz-DTTCompatibility-class),
13
- software_ref
(zzz-DTTCompatibility-class),
13
- software_ref,Software_assnref-method
(zzz-DTTCompatibility-class),
13
- software_ref,Softwares_assnreflist-method
(zzz-DTTCompatibility-class),
13
- Software_ref-class
(zzz-DTTCompatibility-class),
13
- softwareApplication
(zzz-DTTCompatibility-class),
13
- softwareApplication,SoftwareApplications_assnlist-method
(zzz-DTTCompatibility-class),
13
- SoftwareApplication-class
(zzz-DTTCompatibility-class),
13
- softwareApplications_assnlist
(zzz-DTTCompatibility-class),
13
- softwareApplications_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class),
13
- SoftwareApplications_assnlist-class
(zzz-DTTCompatibility-class),
13
- softwareManufacturers_assnreflist
(zzz-DTTCompatibility-class),
13
- softwareManufacturers_assnreflist,Software-method
(zzz-DTTCompatibility-class),
13
- SoftwareManufacturers_assnreflist-class
(zzz-DTTCompatibility-class),
13
- softwares_assnreflist
(zzz-DTTCompatibility-class),
13
- softwares_assnreflist,Hardware-method
(zzz-DTTCompatibility-class),
13

- softwares_assnreflist,Protocol-method
(zzz-DTTCompatibility-class),
13
- softwares_assnreflist,Software-method
(zzz-DTTCompatibility-class),
13
- Softwares_assnreflist-class
(zzz-DTTCompatibility-class),
13
- sourceBioAssays_assnreflist
(zzz-DTTCompatibility-class),
13
- sourceBioAssays_assnreflist,BioAssayMap-method
(zzz-DTTCompatibility-class),
13
- SourceBioAssays_assnreflist-class
(zzz-DTTCompatibility-class),
13
- sourceBioMaterialMeasurements_assnlist
(zzz-DTTCompatibility-class),
13
- sourceBioMaterialMeasurements_assnlist,BioAssayCreation-method
(zzz-DTTCompatibility-class),
13
- sourceBioMaterialMeasurements_assnlist,Hybridization-method
(zzz-DTTCompatibility-class),
13
- sourceBioMaterialMeasurements_assnlist,Treatment-method
(zzz-DTTCompatibility-class),
13
- SourceBioMaterialMeasurements_assnlist-class
(zzz-DTTCompatibility-class),
13
- sourceContact_assnreflist
(zzz-DTTCompatibility-class),
13
- sourceContact_assnreflist,BioSource-method
(zzz-DTTCompatibility-class),
13
- SourceContact_assnreflist-class
(zzz-DTTCompatibility-class),
13
- sourcesQuantitationType_assnreflist
(zzz-DTTCompatibility-class),
13
- sourcesQuantitationType_assnreflist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
13
- SourcesQuantitationType_assnreflist-class
(zzz-DTTCompatibility-class),
13
- spacingBetweenCols
(zzz-DTTCompatibility-class),
13
- spacingBetweenCols,ZoneLayout-method
(zzz-DTTCompatibility-class),
13
- spacingBetweenRows
(zzz-DTTCompatibility-class),
13
- spacingBetweenRows,ZoneLayout-method
(zzz-DTTCompatibility-class),
13
- spacingsBetweenZonesX
(zzz-DTTCompatibility-class),
13
- spacingsBetweenZonesX,ZoneGroup-method
(zzz-DTTCompatibility-class),
13
- spacingsBetweenZonesY
(zzz-DTTCompatibility-class),
13
- spacingsBetweenZonesY,ZoneGroup-method
(zzz-DTTCompatibility-class),
13
- specializedQuantitationType
(zzz-DTTCompatibility-class),
13
- specializedQuantitationType,QuantitationType_assnreflist-method
(zzz-DTTCompatibility-class),
13
- SpecializedQuantitationType-class
(zzz-DTTCompatibility-class),
13
- specializedQuantitationType_ref
(zzz-DTTCompatibility-class),
13
- specializedQuantitationType_ref,QuantitationType_assnreflist-method
(zzz-DTTCompatibility-class),
13
- specializedQuantitationType_ref,QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class),
13
- specializedQuantitationType_ref,SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class),
13
- specializedQuantitationType_ref,TargetQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class),
13
- SpecializedQuantitationType_ref-class
(zzz-DTTCompatibility-class),
13
- species_assn
(zzz-DTTCompatibility-class),

- 13
- species_assn, BioSequence-method
(zzz-DTTCompatibility-class),
13
- species_assn, CompositeGroup-method
(zzz-DTTCompatibility-class),
13
- species_assn, FeatureGroup-method
(zzz-DTTCompatibility-class),
13
- species_assn, ReporterGroup-method
(zzz-DTTCompatibility-class),
13
- Species_assn-class
(zzz-DTTCompatibility-class),
13
- StandardQuantitationType_ref-class
(zzz-DTTCompatibility-class),
13
- startCoord
(zzz-DTTCompatibility-class),
13
- startCoord, MismatchInformation-method
(zzz-DTTCompatibility-class),
13
- strandType
(zzz-DTTCompatibility-class),
13
- strandType, SeqFeatureLocation-method
(zzz-DTTCompatibility-class),
13
- subregions_assnlist
(zzz-DTTCompatibility-class),
13
- subregions_assnlist, SeqFeatureLocation-method
(zzz-DTTCompatibility-class),
13
- Subregions_assnlist-class
(zzz-DTTCompatibility-class),
13
- substrateType_assn
(zzz-DTTCompatibility-class),
13
- substrateType_assn, ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- SubstrateType_assn-class
(zzz-DTTCompatibility-class),
13
- summaryStatistics_assnlist
(zzz-DTTCompatibility-class),
13
- summaryStatistics_assnlist, DerivedBioAssayData-method
(zzz-DTTCompatibility-class),
13
- summaryStatistics_assnlist, MeasuredBioAssayData-method
(zzz-DTTCompatibility-class),
13
- SummaryStatistics_assnlist-class
(zzz-DTTCompatibility-class),
13
- surfaceType_assn
(zzz-DTTCompatibility-class),
13
- surfaceType_assn, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- SurfaceType_assn-class
(zzz-DTTCompatibility-class),
13
- target_assnref
(zzz-DTTCompatibility-class),
13
- target_assnref, BioAssayTreatment-method
(zzz-DTTCompatibility-class),
13
- target_assnref, ImageAcquisition-method
(zzz-DTTCompatibility-class),
13
- Target_assnref-class
(zzz-DTTCompatibility-class),
13
- targetQuantitationType_assnref
(zzz-DTTCompatibility-class),
13
- targetQuantitationType_assnref, Error-method
(zzz-DTTCompatibility-class),
13
- targetQuantitationType_assnref, ExpectedValue-method
(zzz-DTTCompatibility-class),
13
- targetQuantitationType_assnref, PValue-method
(zzz-DTTCompatibility-class),
13
- targetQuantitationType_assnref, QuantitationTypeMap-method
(zzz-DTTCompatibility-class),
13
- TargetQuantitationType_assnref-class
(zzz-DTTCompatibility-class),
13
- technologyType_assn
(zzz-DTTCompatibility-class),
13
- technologyType_assn, FeatureGroup-method

- (zzz-DTTCompatibility-class),
13
- TechnologyType_assn-class
(zzz-DTTCompatibility-class),
13
- temperatureUnit
(zzz-DTTCompatibility-class),
13
- temperatureUnit,Unit_assn-method
(zzz-DTTCompatibility-class),
13
- TemperatureUnit-class
(zzz-DTTCompatibility-class),
13
- TemplateFile-class
(zzz-ARRCompatibility-class),
12
- timeUnit (zzz-DTTCompatibility-class),
13
- timeUnit,Unit_assn-method
(zzz-DTTCompatibility-class),
13
- TimeUnit-class
(zzz-DTTCompatibility-class),
13
- tolerance (zzz-DTTCompatibility-class),
13
- tolerance,ArrayManufacture-method
(zzz-DTTCompatibility-class),
13
- tollFreePhone
(zzz-DTTCompatibility-class),
13
- tollFreePhone,Organization-method
(zzz-DTTCompatibility-class),
13
- tollFreePhone,Person-method
(zzz-DTTCompatibility-class),
13
- topLevelBioAssays_assnreflist
(zzz-DTTCompatibility-class),
13
- topLevelBioAssays_assnreflist,ExperimentDesign-method(zzz-ARRCompatibility-class),
(zzz-DTTCompatibility-class),
13
- TopLevelBioAssays_assnreflist-class
(zzz-DTTCompatibility-class),
13
- transformation
(zzz-DTTCompatibility-class),
13
- transformation,ProducerTransformation_assn-method
(zzz-DTTCompatibility-class),
13
- Transformation-class
(zzz-DTTCompatibility-class),
13
- treatment (zzz-DTTCompatibility-class),
13
- treatment,Treatments_assnlist-method
(zzz-DTTCompatibility-class),
13
- Treatment-class
(zzz-DTTCompatibility-class),
13
- treatments_assnlist
(zzz-DTTCompatibility-class),
13
- treatments_assnlist,BioSample-method
(zzz-DTTCompatibility-class),
13
- treatments_assnlist,BioSource-method
(zzz-DTTCompatibility-class),
13
- treatments_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class),
13
- Treatments_assnlist-class
(zzz-DTTCompatibility-class),
13
- type (zzz-DTTCompatibility-class), 13
- type,ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- type,Measurement-method
(zzz-DTTCompatibility-class),
13
- type,NameValueType-method
(zzz-DTTCompatibility-class),
13
- type,PhysicalArray-method
(zzz-ARRCompatibility-class),
12
- type,TemplateFile-method
(zzz-ARRCompatibility-class),
12
- type,UserAttribute-method
(zzz-ARRCompatibility-class),
12
- type_assn (zzz-DTTCompatibility-class),
13
- type_assn,BioSample-method
(zzz-DTTCompatibility-class),
13

- type_assn, BioSequence-method
(zzz-DTTCompatibility-class),
[13](#)
- type_assn, DatabaseEntry-method
(zzz-DTTCompatibility-class),
[13](#)
- type_assn, DerivedBioAssay-method
(zzz-DTTCompatibility-class),
[13](#)
- type_assn, Hardware-method
(zzz-DTTCompatibility-class),
[13](#)
- type_assn, NodeValue-method
(zzz-DTTCompatibility-class),
[13](#)
- type_assn, Protocol-method
(zzz-DTTCompatibility-class),
[13](#)
- type_assn, Software-method
(zzz-DTTCompatibility-class),
[13](#)
- Type_assn-class
(zzz-DTTCompatibility-class),
[13](#)
- types_assnlist
(zzz-DTTCompatibility-class),
[13](#)
- types_assnlist, CompositeGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- types_assnlist, ExperimentDesign-method
(zzz-DTTCompatibility-class),
[13](#)
- types_assnlist, FeatureGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- types_assnlist, ReporterGroup-method
(zzz-DTTCompatibility-class),
[13](#)
- Types_assnlist-class
(zzz-DTTCompatibility-class),
[13](#)
- unit_assn (zzz-DTTCompatibility-class),
[13](#)
- unit_assn, Measurement-method
(zzz-DTTCompatibility-class),
[13](#)
- Unit_assn-class
(zzz-DTTCompatibility-class),
[13](#)
- unitName (zzz-DTTCompatibility-class),
[13](#)
- unitName, ConcentrationUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitName, DistanceUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitName, MassUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitName, QuantityUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitName, TemperatureUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitName, TimeUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitName, VolumeUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV, ConcentrationUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV, DistanceUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV, MassUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV, QuantityUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV, TemperatureUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV, TimeUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- unitNameCV, VolumeUnit-method
(zzz-DTTCompatibility-class),
[13](#)
- upperLeftX
(zzz-DTTCompatibility-class),
[13](#)
- upperLeftX, Zone-method
(zzz-DTTCompatibility-class),
[13](#)
- upperLeftY

- (zzz-DTTCompatibility-class),
13
- upperLeftY, Zone-method
(zzz-DTTCompatibility-class),
13
- uri (zzz-DTTCompatibility-class), 13
- uri, BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- uri, Database-method
(zzz-DTTCompatibility-class),
13
- uri, DatabaseEntry-method
(zzz-DTTCompatibility-class),
13
- uri, Description-method
(zzz-DTTCompatibility-class),
13
- uri, Hardware-method
(zzz-DTTCompatibility-class),
13
- uri, Image-method
(zzz-DTTCompatibility-class),
13
- uri, Organization-method
(zzz-DTTCompatibility-class),
13
- uri, Person-method
(zzz-DTTCompatibility-class),
13
- uri, Protocol-method
(zzz-DTTCompatibility-class),
13
- uri, Software-method
(zzz-DTTCompatibility-class),
13
- user (accessors), 2
- user, NetAffxResource-method
(NetAffxResource-class), 5
- userAttribute
(zzz-ARRCompatibility-class),
12
- userAttribute, UserAttributes-method
(zzz-ARRCompatibility-class),
12
- UserAttribute-class
(zzz-ARRCompatibility-class),
12
- userAttributes
(zzz-ARRCompatibility-class),
12
- userAttributes, ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- UserAttributes-class
(zzz-ARRCompatibility-class),
12
- userAttributeValue
(zzz-ARRCompatibility-class),
12
- userAttributeValue, UserAttribute-method
(zzz-ARRCompatibility-class),
12
- UserAttributeValue-class
(zzz-ARRCompatibility-class),
12
- value (zzz-ARRCompatibility-class), 12
- value, Control-method
(zzz-ARRCompatibility-class),
12
- value, Datum-method
(zzz-DTTCompatibility-class),
13
- value, Measurement-method
(zzz-DTTCompatibility-class),
13
- value, NameValueType-method
(zzz-DTTCompatibility-class),
13
- value, NodeValue-method
(zzz-DTTCompatibility-class),
13
- value, OntologyEntry-method
(zzz-DTTCompatibility-class),
13
- value, ParameterValue-method
(zzz-DTTCompatibility-class),
13
- value_assn
(zzz-DTTCompatibility-class),
13
- value_assn, FactorValue-method
(zzz-DTTCompatibility-class),
13
- Value_assn-class
(zzz-DTTCompatibility-class),
13
- vector, 4
- version (zzz-ARRCompatibility-class), 12
- version, ArrayDesign-method
(zzz-DTTCompatibility-class),

- 13
- version, ArraySetFile-method
(zzz-ARRCompatibility-class),
12
- version, Database-method
(zzz-DTTCompatibility-class),
13
- version, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13
- version, SoftwareApplication-method
(zzz-DTTCompatibility-class),
13
- version, TemplateFile-method
(zzz-ARRCompatibility-class),
12
- volume (zzz-DTTCompatibility-class), 13
- volume, BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- volumeUnit
(zzz-DTTCompatibility-class),
13
- volumeUnit, Unit_assn-method
(zzz-DTTCompatibility-class),
13
- VolumeUnit-class
(zzz-DTTCompatibility-class),
13
- warningType_assn
(zzz-DTTCompatibility-class),
13
- warningType_assn, Reporter-method
(zzz-DTTCompatibility-class),
13
- WarningType_assn-class
(zzz-DTTCompatibility-class),
13
- width (zzz-DTTCompatibility-class), 13
- width, ArrayGroup-method
(zzz-DTTCompatibility-class),
13
- x (zzz-DTTCompatibility-class), 13
- x, Position-method
(zzz-DTTCompatibility-class),
13
- xclass, 8–10, 11
- xmlValue.XMLAttributeValue, 12
- y (zzz-DTTCompatibility-class), 13
- y, Position-method
(zzz-DTTCompatibility-class),
13
- year (zzz-DTTCompatibility-class), 13
- year, BibliographicReference-method
(zzz-DTTCompatibility-class),
13
- zone (zzz-DTTCompatibility-class), 13
- zone, ZoneLocations_assnlist-method
(zzz-DTTCompatibility-class),
13
- Zone-class
(zzz-DTTCompatibility-class),
13
- zone_assnref
(zzz-DTTCompatibility-class),
13
- zone_assnref, Feature-method
(zzz-DTTCompatibility-class),
13
- zone_assnref, ZoneDefect-method
(zzz-DTTCompatibility-class),
13
- Zone_assnref-class
(zzz-DTTCompatibility-class),
13
- zone_ref (zzz-DTTCompatibility-class),
13
- zone_ref, Zone_assnref-method
(zzz-DTTCompatibility-class),
13
- Zone_ref-class
(zzz-DTTCompatibility-class),
13
- zoneDefect
(zzz-DTTCompatibility-class),
13
- zoneDefect, Adjustments_assnlist-method
(zzz-DTTCompatibility-class),
13
- ZoneDefect-class
(zzz-DTTCompatibility-class),
13
- zoneGroup (zzz-DTTCompatibility-class),
13
- zoneGroup, ZoneGroups_assnlist-method
(zzz-DTTCompatibility-class),
13
- ZoneGroup-class
(zzz-DTTCompatibility-class),
13
- zoneGroups_assnlist

(zzz-DTTCompatibility-class),
13

zoneGroups_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class),
13

ZoneGroups_assnlist-class
(zzz-DTTCompatibility-class),
13

zoneLayout
(zzz-DTTCompatibility-class),
13

zoneLayout,ZoneLayout_assn-method
(zzz-DTTCompatibility-class),
13

ZoneLayout-class
(zzz-DTTCompatibility-class),
13

zoneLayout_assn
(zzz-DTTCompatibility-class),
13

zoneLayout_assn,ZoneGroup-method
(zzz-DTTCompatibility-class),
13

ZoneLayout_assn-class
(zzz-DTTCompatibility-class),
13

zoneLocations_assnlist
(zzz-DTTCompatibility-class),
13

zoneLocations_assnlist,ZoneGroup-method
(zzz-DTTCompatibility-class),
13

ZoneLocations_assnlist-class
(zzz-DTTCompatibility-class),
13

zonesPerX (zzz-DTTCompatibility-class),
13

zonesPerX,ZoneGroup-method
(zzz-DTTCompatibility-class),
13

zonesPerY (zzz-DTTCompatibility-class),
13

zonesPerY,ZoneGroup-method
(zzz-DTTCompatibility-class),
13

zzz-.TypedAssociation-class
(Internal-class), 4

zzz-AffyCompatible-class
(AffyCompatible-class), 3

zzz-ARRCompatibility-class, 12

zzz-DTTCompatibility-class, 13

zzz-NetAffxCompatibility-class, 14